



**KARMAN HEALTHCARE INC.**  
19255 San Jose Ave.  
City Of Industry, CA 91748  
Tel: 626-581-2235 • Fax: 626-581-2335  
Toll Free: 800-805-2762



### 6061 T6 Aircraft-Grade Aluminum

## VIP515 SERIES ORDER FORM

Account Information | Account No.: \_\_\_\_\_ Ship To | \_\_\_\_\_  
Date: \_\_\_\_\_ Name \_\_\_\_\_  
P.O. #: \_\_\_\_\_ Address \_\_\_\_\_  
Buyer: \_\_\_\_\_ City \_\_\_\_\_ St \_\_\_\_\_ Zip \_\_\_\_\_  
Authorized Signature: \_\_\_\_\_ Phone \_\_\_\_\_

**\*Special Note:** HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

### VIP515 Base Chair

MSRP

( ) **VIP515:** Ultra-Lightweight Tilt-In-Space

**\$3,330**

Back **Wheelchair** w/ 20" Flat Free Rear Wheel & Hand Rim, 7"x1" Polyurethane Flat Free Casters, T-6 Aircraft Grade Aluminum 38lb Frame, Swing Away Footrest, Flip-Back Armrests, Breathable Mesh Seat & Back Upholstery, Diamond Black Frame Color, 250lb Weight Capacity, 0°- 35° Tilt Range

( ) **VIP515TP:** Ultra-Lightweight Tilt-In-Space **\$3,150**

Back **Transport Chair** w/ 14" Flat Free Rear Wheel, 7"x1" Polyurethane Flat Free Casters, T-6 Aircraft Grade Aluminum 34lb Frame, Swing Away Footrest, Flip-Back Arms, Breathable Mesh Seat & Back Upholstery, Diamond Black Frame Color, Weight Capacity 250lb, 0°- 35° Tilt Range

❖ **Seat Width:** ( ) 16" ( ) 18"

❖ **Seat Depth:** ( ) 17" Standard

Or ( ) 19" Seat Depth + \$225

### Seat Cushion

- |   |       |
|---|-------|
| ( ) 3" 16x16 Foam                         | \$72  |
| ( ) 3" 18x16 Foam                         | \$72  |
| ( ) 3" 16x16 Foam + Gel                   | \$98  |
| ( ) 3" 18x16 Foam + Gel                   | \$98  |
| ( ) 2 ½" 16x16 Contoured Ergo Memory Foam | \$176 |
| ( ) 2 ½" 18x16 Contoured Ergo Memory Foam | \$176 |

### Armrests

- ( ) Height Adjustable, Flip-back Std

### Hanger Type

- |                                       |            |
|---------------------------------------|------------|
| ( ) Standard Swing Away               | <b>Std</b> |
| ( ) Swing Away Elevating/Articulating | \$278      |

### Optional Equipment

- |   |       |
|---|-------|
| ( ) Substitute Headrest: Adjustable Angle/Length/Height, Rotates 360° | \$238 |
| ( ) Cushion For Detachable Headrest                                   | \$75  |
| ( ) Pelvic Seat Belt w/ composite clamps                              | \$19  |
| ( ) Pelvic Seat Belt, w/ Auto Style Push Button                       | \$24  |
| ( ) AEGIS Anti-Bacterial Seat   | \$56  |
| ( ) AEGIS Anti-Bacterial Back   | \$56  |
| ( ) Anti-Tippers  | \$56  |
| ( ) Heel Loop   | \$18  |

**TOTAL CHAIR PRICE:** \$ \_\_\_\_\_

\*Prices and Specifications Subject to Change without Notice

**COMMENTS:** \_\_\_\_\_