

S-106

Ergonomic Class

ULTRA LIGHTWEIGHT ADJUSTABLE BACKREST

SPECIFICATIONS

S-106-F18

| | |
|---------------------|----------|
| Seat Width | 18 in. |
| Seat Depth | 17 in. |
| Armrest Height | 8 in. |
| Seat Height | 19 in. |
| Back Height | 17 in. |
| Overall Height | 36 in. |
| Overall Open Width | 25 in. |
| Folded Width | 12 in. |
| Overall Length | 39 in. |
| Weight w/o Riggings | 27 in. |
| Weight Capacity | 250 lbs. |
| HCPCS Code | K0004 |

(SUBJECT TO CHANGE WITHOUT NOTICE)



3 EASY ANGLE ADJUSTMENTS FOR THE **IDEAL** POSITION

- Adjustable angle backrest in 3 adjustments - 95°, 100°, 105°
- Fixed armrests and footrest
- High strength frame
- Lightweight (27 lbs.)
- Anti-bacterial, removable/washable upholstery
- Folding precision backrest
- "Tube-in-center" footplate assures side leg support
- 24" Rear polyurethane, high-profile flat-free wheels
- Rear anti-tippers