

VIP-515 Series

World's Only Foldable Lightweight Tilt-in-Space



- Intelligent 0-35 Center-of-Gravity Weight Shifting
- AEGIS Permanently Bonded Antimicrobial Barrier Upholstery Protects from Odor, Staining and Deterioration from Bacteria, Fungus and Other Micro-Organisms
- Bio-Technology : Paint Coating Uses Special Anti-Bacterial Materials
- Exclusive Perfectly Matching Gas Cylinders
- Special Newly Designed Ergonomic Pillows
- Flip-back & Height Adjustable Armrests
- Newly Designed Swing In & Out Footrests
- Weight Capacity: 250 lbs.

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
VIP515TP16	16	17	9 - 11	20	18 + 16	53	24	15	44	34	250
VIP515TP18	18	17	9 - 11	20	18 + 16	53	26	15	44	34	250
VIP515-16	16	17	9 - 11	20	18 + 16	54	26	15	44	38	250
VIP515-18	18	17	9 - 11	20	18 + 16	54	28	15	44	38	250

**Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam Cushion, Ergonomic/Contoured Cushion Reflecting Spoke Guard, Frog Leg Suspension, Elevating Footrests*



Karman Healthcare Inc.
19255 San Jose Ave.
City of Industry, CA 91748
Phone: 626-581-2235 Fax: 626-581-2335
Email: info001@karmanhealthcare.com