



Taking Aluminum to Titanium Heights Titanium Weight at Aluminum price!

ErgoFlight Ultra Lightweight Ergonomic Wheelchair

- 14.5 lb Aluminum Frame. This Manual Wheelchair Competes Against Titanium
- Patented S-Style Ergonomic Seat Frame
- 6061 T-6 Aircraft-grade Aluminum Only 19.8 lbs!
- Companion Brakes Located on Handles for Added Safety and Comfort.
- Built in Silver Aegis Anti-bacterial Seat and Back Cushion
- Folding Backrest/Folding Seat for Transporting in Vehicle or Travel
- Fixed Armrests w/Ergonomic Arm Pads
- Removable Swing-in & Out Footrests
- 20" Rear Polyurethane, High Tread, Flat Free Wheels
- 6" x 1" Polyurethane Front Casters
- 18" x 17", 16" x 17" Seat Width available
- Available Quick Release Rear Wheels
- Frame Color: Pearl Silver
- Weight Capacity: 220 lbs
- Limited Lifetime Warranty on Frame

U.K. ISO 7176/19 Dynamic Crash Test Passed

	seat	seat	armrest	seat	back	overall	overall	folded	overall	weight w/o	Wt.
	width	depth	height	height	height	height	open width	width	length	riggings	Cap.
Model No.	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(lb.)	(lb.)
S-2512F18SS	18	17	8	18.5	17.5	36	24 ¾	15	38	19	220
S-2512F16SS	16	17	8	18.5	17.5	36	22 3⁄4	15	38	19	220
S-2512Q18SS	18	17	8	18.5	17.5	36	24 ¾	15	38	19	220
S-2512Q16SS	16	17	8	18.5	17.5	36	22 3⁄4	15	38	19	220

*Available Options ~ Quick Release wheels, Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam cushions, UV Resistant Canopy, Rigidifying Push Bar, Foldable/Adjustable Headrest

Karman Healthcare Inc.

19255 San Jose Avenue, City of Industry, CA 91748-1418 Order.1-800-805-2762 / Tel.626-581-2235 / Fax.626-581-2335