



R-4800 Aluminum Rollator

Extra Wide Aluminum 4-wheeled Rollators w/8" wheels

Folding/Lightweight aluminum Rollator w/powder coat finish in burgundy or blue

Adjustable height, Hand-grips with locking loop brakes on (4)8" casters

Padded large seat and curved/folding backrest for added comfort + easy to fold and store

Easy-remove under seat storage basket with handle

Easier to maneuver outdoors and indoors than standard Rollators with 8" casters

Weight w/ Basket is 20 lbs

Weight Capacity is 400 lbs

	Seat Width (in.)	seat pad depth (in.)	Seat pad width (in.)	seat height (in.)	backrest height (in.)	overall height (in.)	overall open width (in.)	weight (lb.)	weight capacity (lb.)
R-4800	22	13	18	22	10.5	35~39	29	20	400

** Price & specifications are subject to change without notice.



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335

Filename: R-4700W
Directory: C:\Users\David\Documents
Template: C:\Users\David\AppData\Roaming\Microsoft\Templates\Normal.dotm
Title: E-800 Series
Subject:
Author: Uriel
Keywords:
Comments:
Creation Date: 6/19/2012 3:53:00 PM
Change Number: 5
Last Saved On: 8/2/2012 3:15:00 PM
Last Saved By: David
Total Editing Time: 21 Minutes
Last Printed On: 8/2/2012 3:15:00 PM
As of Last Complete Printing
Number of Pages: 1
Number of Words: 147 (approx.)
Number of Characters: 839 (approx.)