

## Features for S-Ergonomic Series

- Ergonomic, ultra lightweight specialized frame
- Pressure relief, stabilization, and anti-slippage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider/Concaved armpad, providing comfort & support



## S-Ergo 115

Ultra Lightweight Ergonomic Manual Wheelchair

- Medicare code K0004 & K0005
- S-Style Ergonomic Seat
- High strength, weighs only 25 lbs. (w/o footrests)
- Silver 1/4" Anti-Bacterial cushion, removable for washing
- Ergonomic handrims
- Folding backrest/folding seat for easy travel
- Fixed armrest w/ wider concaved armpads
- Swing In & Away footrest allowing maximum safety while entering or exiting the chair
- "Tube-in Center" foot-plate, assures better side leg support
- Pouch for carrying small items attached to upholstery located behind backrest
- 24' Rear spooked Wheel w/polyurethane, high tread, flat free wheels
- 7x1" Polyurethane front casters
- 16"x17", 18"x17", 20"x 17"(250 lbs capacity/26 lbs)
- Frame Color: Pearl Silver & Rose Red
- Upholstery: Black breathable mesh bottom, AEGIS top
- Weight Capacity: 250 lbs

	seat	seat	armrest	seat	back	overall	overall	folded	overall	weight w/o	weight
	width	depth	height	height	height	height	open width	width	length	riggings	capacity
Model No.	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(lb.)	(lb.)
S-115F16	16	17	8	19	17	36	23	12	39	25	250
S-115F18	18	17	8	19	17	36	25	12	39	25	250
S-115F20	20	17	8	19	17	36	27	12	39	25	250

\*Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam Cushion, UV Resistant Canopy, Rigidifying Push-Bar, Reflecting Spoke Guard, Foldable/Adjustable Headrest, Frog Leg Suspension, Magnesium Wheel, Vulcan Grip Wheel, Mag Wheel(Black Plastic)

\*S-115Q – Quick-release rear wheel model, same specifications as above.

