

## ErgoFlight Ergonomic Class

### ULTRA LIGHTWEIGHT AIRCRAFT GRADE T6 ALUMINUM

SPECIFICATIONS	FLIGHT-16	FLIGHT-18
Seat Width	16 in.	18 in.
Seat Depth	17 in.	17 in.
Armrest Height	8 in.	8 in.
Seat Height	18½ in.	18½ in.
Back Height	17½ in.	17½ in.
Overall Width	36 in.	36 in.
Overall Open Width	25 in.	27 in.
Folded Width	15 in.	15 in.
Overall Length	38 in.	38 in.
Weight w/o Riggings	19 lbs.	19 lbs.
Weight Capacity	220 lbs.	220 lbs.
HCPCS Code	K0004	K0004

(SUBJECT TO CHANGE WITHOUT NOTICE)



## GEOMETRY AND TECHNOLOGY TO ITS **BEST...**

- S-Style Ergonomic Seat
- 6061 T-6 Aircraft-grade Aluminum Only 19.8 lbs. (w/o footrests)
- Folding Backrest/ Folding Seat for Transporting in Vehicle or Travel
- Fixed Armrests w/ Concaved Armpads
- 20" Rear Polyurethane, High Tread, Flat Free Wheels
- 6061 T-6 Aircraft-grade Aluminum Only 14.5 lb
- Patented S-Style Ergonomic Seat Frame
- Built in Silver Aegis Anti-bacterial Cushion
- 18" x 17" Seat Width
- Frame Color: Pearl Silver