

## LT-2000 Transport Class

# TRANSPORTER

## DETACHABLE ARMREST

SPECIFICATIONS	LT-2017	LT-2019
Seat Width	17 in.	19 in.
Seat Depth	16 in.	16 in.
Seat Height	19½ in.	19½ in.
Back Height	19 in.	19 in.
Overall Height	40 in.	40 in.
Overall Open Width	21¼ in.	23¼ in.
Weight w/o Riggings	19 lbs.	19 lbs.
Weight Capacity	250 lbs.	250 lbs.
HCPCS Code	E1038	E1038

(SUBJECT TO CHANGE WITHOUT NOTICE)



## FIND THE **PERFECT** FIT FOR YOU...

- Padded fixed armrests
- Foldable w/ folding backrest
- 17" or 19" Seat Width of nylon
- Weight of 19 lbs. w/capacity of 250 lbs.
- Colors: Black, Blue, or Burgundy
- High Strength – Aluminum Frame