

Karman Healthcare

PRODUCT CATALOG 2013

Unlimited Mobility Solutions for Active and Standard Users

Solutions for Active and Standard Users

Healthcare Inc. All Rights Reserved.

Table of Contents

Terms and Conditions - 2

Ergonomic Wheelchairs - 3

Ergo Flight - 4

S-Ergo Lite - 5

S-Ergo 105 - 6

S-Ergo 106 - 7

S-Ergo 115 - 8

S-Ergo ATX - 9

S-Ergo 125 - 10

S-Ergo 305 - 11

Lightweight Wheelchairs - 12

LT-700T - **13**

LT-800T - 14

802-DY - 15

KM-3520 - 16

KM-802F-SOMA - 17

KM-9020L - 18

LT-K5 - **19**

LT-980 - 20

Standard Wheelchairs - 21

KN-700T - 22

KN-800T - 23

Bariatric Wheelchairs - 24

KM-BT10 - 25

KM-8520 - 26

KN-920W - 27

KN-922W - 28

KN-924W, KN-926W, KN-928 - 29

T-920/ T-922 - 30

KN-880-WE - 31

KM-5000F - 32

MVP-502-MS - 33

Transport Wheelchairs - 34

T-2700 - 35

LT-1000 - 36

LT-2000 - 37

S-115-TP - 38

S-Ergo Lite - 39

T-2000 - **40**

T-900 - **41**

1-900 - 41

KM-TV10B - 42

MVP-502 - **43**

VIP-515 - 44

KM-5000 - **45**

Stand Up Wheelchairs - 46

XO-101 - **47**

XO-202 - 48

Active Wheelchairs - 49

AT-60 - **50**

S-Ergo ATX - 51

Rollators - 52

R-3600 / R-4100 - 53

R-4100 / R-4200 - 54

R-4600 / R-4602 - 55

R-4608 / R-4700 - 56

R-4800 / KW-100 - 57

Daily Aids - Accessories - 58



Ergonomic Wheelchairs

- Ergo Flight 4
- S-Ergo Lite 5
- S-Ergo 105 **6**
- S-Ergo 106 **7**
- S-Ergo 115 8
- S-Ergo ATX 9
- S-Ergo 125 **10**
- S-Ergo 305 **11**

Ergo Flight Wheelchair



Product Features

One of the lightest Ergonomic wheelchairs on the market, weighing in at only 19.8 lbs.

The Ergo Flight combines a lightweight 6061 T6 aluminum frame for ease of use with a patented S-shaped seating frame to provide better comfort and more security when in motion.

Color: Pearl Silver

- User operated braking system
- S-Shaped seat
- AEGIS anti-bacterial seat cover
- Ultra light T6 alumininum frame
- Quick Release Available

Product Specifications	Ergo Flight-16	Ergo Flight-18	Ergo Flight-16Q	Ergo Flight-18Q
Seat Width:	16 in.	18 in.	16 in.	18 in.
Seat Depth:	17 in.	17 in.	17 in.	17 in.
Seat Height:	18.5 in.	18.5 in.	18.5 in.	18.5 in.
Back Height:				
Overall Height:	36 in.	36 in.	36 in.	36 in.
Overall Folded Width:	15 in.	15 in.	15 in.	15 in.
Overall Length:	38 in.	38 in.	38 in.	38 in.
Weight:	19.5 lbs.	19.5 lbs.	19.5 lbs.	19.5 lbs.
Weight Capacity:	220 lbs.	220 lbs.	220 lbs.	220 lbs.
HCPCS Code:	N/A	N/A	N/A	N/A

Ergo Lite Wheelchair



Product Features

What do you get when you combine an ultra light 19 lb aluminum frame and the comfort of Ergonomic seating?

You get one of our best selling transport chairs, the Ergo Lite! This transporter is perfect for those who need quality comfort at an affordable price.

Color: Pearl Silver

- Lightest transporter on the market.
- S-Shape seating provides weight distribution for comfort.
- User controlled parking brakes on handle grips.
- Folding backrest for easy transportation.

Product Specifications	Ergo Lite-16	Ergo Lite-18
Seat Width:	16 in.	18 in.
Seat Depth:	17 in.	17 in.
Seat Height:	19 in.	19 in.
Back Height:		
Overall Height:	36 in.	36 in.
Overall Folded Width:	12 in.	12 in.
Overall Length:	39 in.	39 in.
Weight:	18 lbs.	18 lbs.
Weight Capacity:	220 lbs.	220 lbs.
HCPCS Code:	N/A	N/A

S-Ergo 105 Wheelchair



Product Features

High strength standard wheelchair in a simple lightweight package, with all the comfort of the Ergo series.

Color: Pearl Silver Color: Rose Red

- S-Shaped Ergo seat provides even distribution of weight for extra comfort
- Fixed footrest
- Extending push to lock manual brakes
- Upholstery: Silver Aegis over Black Breathable Mesh bottom

Product Specifications	S-105-F16	S-105-F18
Seat Width:	16 in.	18 in.
Seat Depth:	17 in.	17 in.
Seat Height:	19 in.	19 in.
Back Height:		
Overall Height:	36 in.	36 in.
Overall Folded Width:	23 in.	23 in.
Overall Length:	39 in.	39 in.
Weight:	27 lbs.	27 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	K0004	K0004

S-Ergo 106 Wheelchair



Product Features

One step up from the Ergo 105 wheelchair. Adds a reclining function that can be operated while the user is in the chair.

Comfort never felt so standard!

- Same light weight 6061 T6 aluminum frame that is used in the Ergo-105
- "Tube-in Center" footplate, assures better side leg support
- Upholstery: Silver ¼" Anti-Bacterial cushion AEGIS, removable, machine washable.

roduct Specifications	S-106-F18
Seat Width:	18 in.
Seat Depth:	17 in.
Seat Height:	19 in.
Back Height:	
Overall Height:	36 in.
Overall Folded Width:	12 in.
Overall Length:	39 in.
Weight:	27 lbs.
Weight Capacity:	250 lbs.
HCPCS Code:	K0004

S-Ergo 115 Wheelchair



Product Features

Light weight Ergonomic wheelchair with folding backrest and a removable axil for easy storage. Perfect for users where storage space is limited.


All the comfort of the Ergo series S-Shaped seats.

Color: Pearl Silver Color: Rose Red

- Swing in and away footrests
- High Strength frame, weighing only 25 lbs.
- 7x1" Polyurethane front casters for a smoother ride
- Silver 1/4" Aegis Anti-Bacterial Upholstery,
- Ergonomic Handrims & S-Style Ergonomic Seat

Product Specifications	S-115F16	S-115F18	S-115F20
Seat Width:	16 in.	18 in.	20 in.
Seat Depth:	17 in.	17 in.	17 in.
Seat Height:	19 in.	19 in.	19 in.
Back Height:			
Overall Height:	36 in.	36 in.	36 in.
Overall Folded Width:	12 in.	12 in.	12 in.
Overall Length:	38 in.	38 in.	38 in.
Weight:	25 lbs.	25 lbs.	25 lbs.
Weight Capacity:	250 lbs.	250 lbs.	250 lbs.
HCPCS Code:	K0004	K0004	K0004

S-Ergo ATX Wheelchair



Product Features

The Ferrari of the Ergo series! Made for those who need a chair to keep up with an active lifestyle. The S-Ergo ATX transfers your energy to the street with ultimate style, comfort, and expression.

*Optional high performance Spinergy wheels available.

Color: Lamborghini Verde Ithaca

- Light weight aluminum frame.
- S-Shaped seat for keeping your body secure while moving.
- Built for active use!

Product Specifications	S-ERGO-ATX 14/16/18	
Seat Width:	18 in.	
Seat Depth:	18 in.	
Seat Height:	14 in.	
Back Height:		
Overall Height:	var	
Overall Folded Width:	var	
Overall Length:	var	
Weight:	15.4 lbs	
Weight Capacity:	275 lbs.	
HCPCS Code:	N/A	

S-Ergo 125 Wheelchair



Product Features

Ultra lightweight wheelchair with the S-115 model standards, plus flip back armrests and color options.

Great for those who may have trouble getting in and out of their chair. Comes with all the bells and whistles of the Ergo series.

Color: Pearl Silver Color: Rose Red

- Swing-away detachable footrests with heel loop
- High-grade/Padded nylon flame retardent upholstery
- Push to lock wheel breaks
- Fixed padded full armrests with alloy side guard
- 7" Flat free polyurethane casters

Product Specifications	S-125F16	S-125F18
Seat Width:	16 in.	18 in.
Seat Depth:	17 in.	17 in.
Seat Height:	18.5 in.	18.5 in.
Back Height:		
Overall Height:	36 in.	36 in.
Overall Folded Width:	15 in.	15 in.
Overall Length:	38 in.	38 in.
Weight:	19.5 lbs.	19.5 lbs.
Weight Capacity:	220 lbs.	220 lbs.
HCPCS Code:	N/A	N/A

S-Ergo 305 Wheelchair



Product Features

Durability and ease of use never felt so perfect together. The S-Ergo 305 combines the durability of standard wheelchairs with the comfort of the Ergonomic series.

*2 optional colors avalible.

- Folds to a compact 12 inches!
- Optional Cushion colors available.

roduct Specifications	S-ERGO-305
Seat Width:	16 in.
Seat Depth:	17 in.
Seat Height:	19 in.
Back Height:	
Overall Height:	38 in.
Overall Folded Width:	12 in.
Overall Length:	43 in.
Weight:	29 lbs.
Weight Capacity:	250 lbs.
HCPCS Code:	K0005



Lightweight Wheelchairs

LT-700T **13**

LT-800T **14**

802-DY **15**

KM-3520 **16**

KM-802F-SOMA 17

KM-9020L **18**

LT-K5 **19**

LT-980 **20**

LT-700T Wheelchair



Product Features


The Karman Standard Lightweight Wheelchair Lt 700 is economical with detachable armrests and swing-away footrests. It can be customized with pneumatic tires and elevating footrests.

Color: Black

- Detachable / Padded desk length armrests with allow side guard
- Swing away/ Length adjustable footrests with heel loop footplates
- Push to lock wheel brakes 24 x 1"
- Mag wheels with polyurethane tires
- Extra-durable forks 8 x 1 flat free front casters

Product Specifications	LT-700NT	LT-700T
Seat Width:	16 in.	18 in.
Seat Depth:	16 in.	16 in.
Armrest Height		
Seat Height:	17/19 in.	17/19 in.
Back Height:		
Overall Width:	36 in.	36 in.
Overall Folded Width:	11 in.	11 in.
Overall Length:	31 in.	31 in.
Weight:	36 lbs.	36 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	K0003	K0003

LT-800T Wheelchair



Product Features

Standard lightweight wheelchair with detachable arm rests. Built to last, with a reinforced steel frame and alloy side guards.

Color: Silver Vein

- Extra durable forks
- Swing-away footrests
- High-grade nylon flameproof seat.
- Push to lock wheel brakes
- Fixed armrests with alloy guards

Product Specifications	LT-800NT	LT-800T
Seat Width:	16 in.	18 in.
Seat Depth:	16 in.	16 in.
Armrest Height		
Seat Height:	19 in.	19 in.
Back Height:	16 in.	16 in.
Overall Width:	36 in.	36 in.
Overall Folded Width:	11 in.	11 in.
Overall Length:	31 in.	31 in.
Weight:	34 lbs.	34 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	K0003	K0003

802-DY Wheelchair



Product Features

This Karman 802-DY Lightweight Wheelchair is a folding high strength aluminum frame with a black powder coat finish. It features comfortable padded flip-back, desk length armrests with plastic side panels.

Swing-away detachable footrests with heel strap and adjustable length anti-slip footplates.

- Folding high strength aluminum frame
- Padded flip-back desk armrests
- Swing-away detachable footrests with heel strap
- 24" x 1" mag wheels with polyurethane tires
- Push to lock wheel brakes

Product Specifications	802-DY-16	802-DY-18
Seat Width:	16 in.	18 in.
Seat Depth:	16 in.	16 in.
Armrest Height		
Seat Height:	17/19 in.	17/19 in.
Back Height:		
Overall Width:	36 in.	36 in.
Overall Folded Width:	11 in.	11 in.
Overall Length:	31 in.	31 in.
Weight:	36 lbs.	36 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	K0003	K0003

KM-3520 Wheelchair



Product Features

This wheelchair features side panels that also become a dinner tray! The wheelchair frame is constructed of 6061 T-6 aircraft-grade aluminum for strength and light weight.

Comes standard with folding backrest for easy storage and removable swingaway footrests with velcro mount leg strap.

Color:Black

- •2-in-1 dinner tray/side panel
- Adjustable table depth
- Folding Backrests
- Padded flip-back armrests

Product Specifications	KM-3520-16	KM-3520-18
Seat Width:	16 in.	18 in.
Seat Depth:	16 in.	16 in.
Seat Height:		
Back Height:	16 in.	16 in.
Overall Width:	38 in.	38 in.
Overall Open Width:	26 in.	28 in.
Overall Height:	38 in.	38 in.
Weight:	29 lbs.	29 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	K0004	K0004

KM-802F Wheelchair







This wheelchair is constructed of 6061 T-6 aircraft-grade aluminum for strength and light weight.

It features flip-back armrests for easy transfers and removable swingaway footrests with velcro mount leg strap. High-grade flame retardant, breathable nylon upholstery.

- 6061 T-6 Aircraft-grade aluminum frame
- Padded full sized/flip-back armrests with side panels
- Swing-away, detachable, height adjustable footrests with leg strap
- Dual axle for adjusting seat height to 18 or 19 in
- 24" x 1" spoke wheels with high profile

Product Specifications	KM-802F16	KM-802F18	KM-802Q16	KM-802Q18
Seat Width:	16 in.	18 in.	16 in.	18 in.
Seat Depth:	17 in.	17 in.	17 in.	17 in.
Armrest Height				
Seat Height:	18 in.	18 in.	18 in.	18 in.
Back Height:	12/18 in.	12/18 in.	12/18 in.	12/18 in.
Overall Height:	38 in.	38 in.	38 in.	38 in.
Overall Open Width:	23 in.	25 in.	23 in.	25 in.
Overall Length:	42 in.	42 in.	42 in.	42 in.
Weight:	30 lbs.	30 lbs.	30 lbs.	30 lbs.
Weight Capacity:	250 lbs.	250 lbs.	250 lbs.	250 lbs.
HCPCS Code:	K0004	K0004	K0004	K0004

KM-9020L Wheelchair



Product Features

The Karman Ultralightweight Compact KM-9020L Wheelchair is easy to lift and its small size makes it ideal for travel.

Weighing only 25 lbs., the chair is easy to store and transport.

Color: Black Color: Champagne

- Folding backrests
- Fixed armrests
- Fixed footrests with velcro mount leg strap
- •6" flat free front casters

Product Specifications	KM-9020L
Seat Width:	18 in.
Seat Depth:	15 in.
Seat Height:	
Back Height:	19 in.
Overall Open Width:	25 in.
Overall Height:	37 in.
Weight:	25 lbs.
Weight Capacity:	220 lbs.
HCPCS Code:	K0004

LT-K5 Wheelchair



Product Features

This wheelchair is constructed of 6061 T-6 aircraft-grade aluminum for strength and light weight.

It features flip-back armrests for easy transfers and removable swingaway footrests with velcro mount leg strap. High-grade flame retardant, breathable nylon upholstery.

- 6061 T-6 Aircraft-grade aluminum frame
- Padded full sized/flip-back armrests with side panels
- Swing-away, detachable, height adjustable footrests with leg strap
- Dual axle for adjusting seat height to 18 or 19 in
- 24" x 1" spoke wheels with high profile

Product Specifications	LT-K5-16	LT-K5-18
Seat Width:	16 in.	18 in.
Seat Depth:	17 in.	17 in.
Armrest Height		
Seat Height:	18 in.	18 in.
Back Height:	12/18 in.	12/18 in.
Overall Height:	38 in.	38 in.
Overall Open Width:	23 in.	25 in.
Overall Length:	42 in.	42 in.
Weight:	30 lbs.	30 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	K0004	K0004

LT-980 Wheelchair



Product Features


Wheelchair is easy to lift and its small size makes it ideal for travel. Weighing only 24 lbs., the chair is easy to store and transport.

The chair folds down the middle like a standard wheelchair and the backrest also folds down to make the wheelchair even more compact.

Color: Black

- Desk-length non-removable armrests
- Height adjustable footplate
- Folding backrest
- Footrests with velcro mount leg strap

Product Specifications	LT-980	LT-980-BD
Seat Width:	18 in.	18 in.
Seat Depth:	15 in.	15 in.
Armrest Height:		
Seat Height:	19 in.	19 in.
Back Height:	19 in.	19 in.
Overall Open Width:	26 in.	26 in.
Overall Height:	41 in.	41 in.
Weight:	24 lbs.	24 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	K0004	K0004



Standard Wheelchairs

KN-700T **13**

KN-800T **14**

KN-700T Wheelchair



Product Features

This wheelchair provides comfort and durability with ease and comes with high grade retardant upholstery.

High quality components allow this wheelchair to thrive in any environment and/or terrain. Detachable armrests and swing away footrests top off the list of great quality components included.

Color: Black

- Detachable padded armrests
- Seat height adjustable
- Push to lock wheel brakes
- High grade retardant upholstery
- Swing away footrests

Product Specifications	KN-700-NT	KN-700T
Seat Width:	16 in.	18 in.
Seat Depth:	16 in.	16 in.
Seat Height:	17/19 in.	17/19 in.
Overall Height:	24 in.	24 in.
Overall Open Width:	24 in.	26 in.
Overall Length:	31 in.	31 in.
Weight:	39 lbs.	39 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	K0001	K0001

KN-800T Wheelchair



Product Features

The Karman essential wheelchair is a standard wheelchair with padded fixed armrests and swingaway, detachable footrests.

It comes equipped with adjustable toggle wheel locks, the length of the footplates is also adjustable.

Color: Black

- Folding powder coating steel frame
- Padded fixed armrests
- Swingaway detachable footrests
- Adjustable toggle wheel locks
- Adjustable length footplates

Product Specifications	KN-800-NT	KN-800T
Seat Width:	16 in.	18 in.
Seat Depth:	16 in.	16 in.
Seat Height:	19 in.	19 in.
Back Height:		
Overall Width:	36 in.	36 in.
Overall Open Width:	23 in.	25 in.
Overall Length:	31 in.	31 in.
Weight:	37 lbs.	37 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	K0001	K0001



Bariatric Wheelchairs

- KM-BT10 25
- KM-8520 **26**
- KN-920W **27**
- KN-922W **28**
- KN-924W, KN-926W, KN-928W **29**
 - T-920/T-922 **30**
 - KN-880-WE **31**
 - MVP-502-MS **32**
 - KM-5000F 33

KM-BT10 Wheelchair



Product Features

BT10 Wheelchair includes a wide range of adjustability features and seat sizes.

The adjustable axle plate allows for center-of-gravity adjustment for optimal performance while the backrest has five angle adjustments available.

Color: Diamond Black

- Unique swing away anti tipper design
- Foldable and robust frame design
- Energy saving self propel mechanism
- Weight supports up to 550 lbs.

Product Specifications	KM-BT10-22W-24W	KM-BT10-26W	KM-BT10-28W-30W
Seat Width:	22 in. / 24 in.	22 in. / 24 in.	28 in. / 30 in.
Seat Depth:	18 / 20 in.	18 - 20 / 22 - 24 in.	22 - 24 in.
Armrest Height:			
Seat Height:	15.5 in 19.5 in.	15.5 in 19.5 in.	15.5 in 19.5 in.
Back Height:	16 - 20 in.	16 - 20 in.	16 - 20 in.
Overall Width:	31 - 33 in.	25 in.	25 in.
Overall Folded Width:	19 in.	15 in.	15 in.
Overall Length:	47 in.	47 in.	47 in.
Weight:	63 - 66 lbs.	63 - 66 lbs.	63 - 66 lbs.
Weight Capacity:	550 lbs.	550 lbs.	550 lbs.
HCPCS Code:	K0007	K0007	K0007

KM-8520W Wheelchair



Product Features

This wheelchair is made of T-6 Aircraft grade aluminum, also includes a folding backrest and swing away detachable footrests with velcro mount leg straps.

More features include heavy duty flat free spoke wheels with polyurethane, made to last you a very long time.

Color: Black

- Swing away detachable footrests
- Flip back armrests
- Adjustable toggle wheel locks
- High grade breathable nylon seat
- Anti slip footplates

Product Specifications	KM08520W20HA	KM08520W20HA-E	KM08520-22W
Seat Width:	20 in.	20 in.	22 in.
Seat Depth:	18 in.	18 in.	18 in.
Color:			
Seat Height:	20 in.	20 in.	20 in.
Back Height:			
Overall Height:	36 in.	36 in.	36 in.
Overall Open Width:	30 in.	30 in.	30 in.
Overall Length:	41 in.	41 in.	41 in.
Weight:	35 lbs.	35 lbs.	35 lbs.
Weight Capacity:	350 lbs.	350 lbs.	350 lbs.
HCPCS Code:	K0007	K0007	K0007

KN-920W Wheelchair



Product Features

A wheelchair needs to be built with durability, portability, and comfort in mind. All these attributes are already synonymous with Karman Healthcare and it's wheelchairs.

The KN-920W Bariatric wheelchair (meaning it has an extra wide seat) fills the gap between the low- and high-end markets, providing an affordable alternative to other offerings.

Color: Chrome

- High quality steel construction
- Detachable padded armrests
- Double padded vinyl upholstery
- Swing away footrests composite footplates

Product Specifications	KN-920W	KN-920WE	
Seat Width:	20 in.	20 in.	
Seat Depth:	16 in.	16 in.	
Seat Height:	18 in.		
Back Height:	18 in.	18 in.	
Overall Width:			
Overall Folded Width:	28 in.	28 in.	
Overall Length:	47 in.	47 in.	
Weight:	41 lbs.	41 lbs.	
Weight Capacity:	300 lbs.	300 lbs.	
HCPCS Code:	K0007	K0007	

KN-922W Wheelchair



Product Features

This wheelchair is the sister model of the T-922W, except this one comes with a wide 22 in. seat and a padded detachable desk length armrests.

High grade flame retardant, breathable nylon upholstery and Dual cross bars.

Color: Black

- Swing away detachable footrests
- Flip back armrests
- Adjustable toggle wheel locks
- High grade breathable nylon seat
- Anti slip footplates

Product Specifications	KN-922W	KN-922WE
Seat Width:	20 in.	20 in.
Seat Depth:	18 in.	18 in.
Color:		
Seat Height:	20 in.	20 in.
Back Height:		
Overall Height:	36 in.	36 in.
Overall Open Width:	30 in.	30 in.
Overall Length:	41 in.	41 in.
Weight:	35 lbs.	35 lbs.
Weight Capacity:	350 lbs.	350 lbs.
HCPCS Code:	K0007	K0007

KN-924,926,928 Wheelchair





STEEL HAND RIMS

Product Features

Built from high quality steel, the KN-W series of wheelchairs are built for durability and function.



Color: Black, Silver Vein

- Steel hand rims for durability
- Detachable padded armrests
- Adjustable seat height
- Double cross braces for support
- Folds compact for storage

Product Specifications	KN-924W	KN-926W	KN-928W
Seat Width:	24 in.	26 in.	28 in.
Seat Depth:	18 in.	20 in.	20 in.
Color:			
Seat Height:	20 in.	19 in.	19 in.
Back Height:		18 in.	
Overall Width:	32 in.	37 in.	37 in.
Weight:	56 lbs.	87 lbs.	89 lbs.
Weight Capacity:	400 lbs.	450 lbs.	600 lbs.
HCPCS Code:	K0007	K0007	K0007

T-920/T-922 Wheelchair





ERGONOMIC ARMRESTS



This Deluxe Bariatric Transport Wheelchair adds to our line of high strength dual cross braces. It supports up to 450 lbs due to its reinforced steel frame. It's more than heavy duty, its peace of mind to know the durability of the wheelchair is there to stand the test of time.



Color: Red

- Dual Cross Bars
- •12.5 in. polyurethane tires
- •450 lb weight capacity
- High-grade flame retardant, breathable nylon upholstery

roduct Specifications	T-920	T-922
Seat Width:	20 in.	22 in.
Seat Depth:	18 in.	18 in.
Seat Height:	21 in.	21 in.
Back Height:		
Overall Width:	25 in.	27 in.
Weight:	40 lbs.	44 lbs.
Weight Capacity:	450 lbs.	450 lbs.
HCPCS Code:	E1038	E1038

KN-880-WE Wheelchair



Product Features

One of the higher end models of our Bariatric chairs. The KN-880 combines the reliability and strength of our steel frame design with the comfort of a Cadillac!

Designed for those who need extra support and comfort. Adds neck and leg support, plus the luxury of a reclining chair.



Color: Chrome

- Fully reclining to 180 degrees
- Detachable elevating footrests
- Adjustable length leg supports
- Folding chrome plated steel frame

roduct Specifications	KN-880NE	KN-880E	KN-880WE
Seat Width:	16 in.	16 in.	20 in.
Seat Depth:	16 in.	16 in.	16 in.
Color:			
Seat Height:	20 in.	20 in.	20 in.
Back Height:	29 in. +	29 in. +	29 in. +
Overall Width:	23 in.	25 in.	27 in.
Weight:	50 lbs.	50 lbs.	52 lbs.
Weight Capacity:	250 lbs.	250 lbs.	300 lbs.
HCPCS Code:	K0003	K0003	K0003

KM-5000F Wheelchair



Product Features

One of only two lightweight recliners available in the USA. Built with a lightweight aluminum frame, this ultra lightweight reclining wheelchair features fully adjustable back rest and leg rests.



Color: Black

- Reclines from 90-163 degrees
- Folding 6061 T-6 aircraft aluminum frame
- Adjustable length leg & foot rests
- Toggle wheel locks

Product Specifications	KM5000-16	KM5000-18	KM5000-20	KM5000-22
Seat Width:	16 in.	18 in.	20 in.	22 in.
Seat Depth:	18 in.	18 in.	18 in.	18 in.
Seat Height:	20 in.	20 in.	20 in.	20 in.
Back Height:	33 in.	33 in.	33 in.	33 in.
Overall Width:	25 in.	27 in.	29 in.	31 in.
Weight:	36 lbs.	36 lbs.	38 lbs.	38 lbs.
Weight Capacity:	250 lbs.	250 lbs.	300 lbs.	350 lbs.
HCPCS Code:	K0004+E1226+E	E0966		

MVP-502-MS Wheelchair



Product Features

Ergonomic V-Seat automatic recliner. The next generation of bariatric chairs. Ultimate comfort with the luxury of an electronic reclining feature. Leg and back support thanks to ergonomic upholstery, which provides pressure relieving comfort throughout the body.



Color: Chrome

- 3D elevating leg rest support
- V-Seating for body pressure relief
- Motorized reclining feature
- Aegis anti-bacterial upholstery

Product Specifications	KN-920W	KN-920WE	KN-920WE
Seat Width:	16 in.	18 in.	20 in.
Seat Depth:	18 in.	18 in.	18 in.
Seat Height:			
Back Height:	19 in.	19 in.	19 in.
Overall Width:		26 in.	28 in.
Overall Folded Width:	14 in.	15 in.	16 in.
Overall Length:	47 in.	47 in.	47 in.
Weight:	36 lbs.	36 lbs.	36 lbs.
Weight Capacity:	250 lbs.	250 lbs.	250 lbs.
HCPCS Code:	K0003/E1226	K0003/E1226	K0003/E1226



Transport Wheelchairs

T-2700 - 35

LT-1000 - 36

LT-2000 - **37**

S-115-TP - 38

Ergo Lite - 39

T-2000 - 40

T-900 - **41**

KM-TV10B - 42

MVP-502-TP - 43

VIP-515-TP - 44

KM-5000-TP - 45

T-2700 Transport Chair



Product Features

Standard transport wheelchair built from a durable steel frame. Foldable with folding backrest for convenient storage.

Comes in both 17 and 19 inch seat sizes.

Color: Black

- · Foldable w/ folding backrest
- · Capacity of 250 lbs!
- Swing-away footrests w/ heel loop
- · Durable steel frame
- · Black finish

Product Specifications	T-2700N	T-2700	
Seat Width:	17 in.	19 in.	
Seat Depth:	16 in.	16 in.	
Color:			
Seat Height:	20 in.	20 in.	
Back Helght:	17 in.		
Overall Width:	25.75 in.	27.75 in.	
Weight:	27 lbs.	27 lbs.	
Weight Capacity:	250 lbs.	250 lbs.	
HCPCS Code:	E1038	E1038	

LT-1000 Transport Chair



Product Features

The LT-1000 is the next series of aluminum high strength, ultra lightweight hand brake transport chairs. Equipped with seat belt to keep user secure. Available in both blue and burgundy color schemes.





Color: Blue

Burgundy

- Ultra lightweight aluminum frame
- Seat belt for user security
- Flame retardant upholstery
- Manual lock brakes
- Only 23 lbs!

Product Specifications	LT-1000-HB	
Seat Width:	19 in.	
Seat Depth:	16 in.	
Overall Open Width:		
Seat Height:	19 in.	
Back Height:	19 in.	
Overall Height:	21.5 in.	
Weight:	23 lbs.	
Weight Capacity:	250 lbs.	
HCPCS Code:	E1038	

LT-2000 Transport Chair



Product Features

The LT-2000 uses an ultra lightweight aluminum frame that folds into a compact form for easy storage, all under 20 lbs! The LT-2000 is the basic version of the LT-1000, for those who don't need the extra bells and whistles like user controlled brakes.



- Padded fixed armrests
- Foldable w/ folding backrest
- Swing away footrests
- 17 in. or 19 in. Seat width of nylon
- High Strength Aluminum Frame

roduct Specifications	LT-2017	LT-2019
Seat Width:	17 in.	19 in.
Seat Depth:	16 in.	16 in.
Seat Height:	19.5 in.	19.5 in.
Back Height:	19 in.	19 in.
Overall Height:	40 in.	40 in.
Overall Open Width:	21.4 in.	23.4 in.
Overall Length:	in.	in.
Weight:	19 lbs.	19 lbs.
Weight Capacity:	250 lbs.	250 lbs
HCPCS Code:	E1038	E1038

S-115 Transport Chair







Scientifically proven to relieve pressure throughout the body, this chair will revolutionize your life! It was designed to highlight the beauty of a customized chair.

Utilizes ergonomics in its seating design accompanied by a breathable anti-bacterial and removable upholstery.





Color: Pearl Silver Color: Rose Red

- S-Style Ergonomic Seat
- Fixed armrests w/ concaved armpads
- Extending push to lock brakes
- Back pouch for extra baggage
- Black breathable mesh bottom

roduct Specifications	S-115TP16	S-115TP18	S-115TP20
Seat Width:	16 in.	18 in.	16 in.
Seat Depth:	17 in.	17 in.	17 in.
Seat Height:	19 in.	19 in.	19 in.
Back Height:	16 in.	16 in.	16 in.
Overall Width:	23 in.	24 in.	25 in.
Weight:	22 lbs.	22 lbs.	23 lbs.
Weight Capacity:	250 lbs.	250 lbs.	250 lbs.
HCPCS Code:	N/A	N/A	N/A

Ergo Lite Transport Chair

3 STAGE BRAKING

FOLDING BACKREST







Product Features

What do you get when you combine an ultra lightweight 19lb aluminum frame with ergonomic seating and cutting edge design? You get the Ergo Lite, one of our best selling lightweight transport chairs.

Lightest transporter on the market!

Color: Pearl Silver

- S-Style ergonomic seat
- Aegis anti-bacterial cushion
- Folding backrest/seat
- Weight capasity of 220 lbs!

roduct Specifications	ErgoLite-16	ErgoLite-18
Seat Width:	16 in.	18 in.
Seat Depth:	17 in.	17 in.
Seat Height:	19 in.	19 in.
Back Height:	17 in.	19 in.
Overall Height:		
Overall Open Width:	23 in.	25 in.
Overall Length:	39 in.	39 in.
Weight:	18 lbs!	18 lbs!
Weight Capacity:	220 lbs.	250 lbs.
HCPCS Code:	N/A	N/A

T-2000 Transport Chair



Product Features

A lightweight compact transport wheelchair with fixed full length armrests & footrest with Velcro mount leg strap.

The frame folds for easy traveling and/or storage and weighs only 23 lbs with the footrest.

Color: Black

- Folding backrest for travel
- Fixed armrest & footrest
- Only 23 lbs. with footrest
- Lightweight steel frame
- Black breathable mesh bottom

Product Specifications	T-2017	T-2019
Seat Width:	17 in.	19 in.
Seat Depth:	16 in.	16 in.
Seat Height:	19 in.	19 in.
Back Height:	19 in.	19 in.
Overall Height:	40 in.	40 in.
Overall Open Width:	22.5 lbs.	23.5 lbs.
Weight:	23 lbs.	23 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	E1038	E1038

T-900 Transport Chair



Product Features

This wheelchair is a deluxe steel bariatric trasport chair that includes dual cross bars for added support. Also features swing away footrests w/ heel loop and manual lock brakes.

High grade flame retardant, breathable nylon upholstery rounds off this trasnport chair that exudes durability and support.

Color: Red

- Padded detachable armrests
- Dual cross bars
- Swing away footrests w/ heel loop
- Manual lock brakes

Seat Width:	20 in.	22 in.
Ceat Depth:	18 in.	18 in.
eat Height:		
Back Height:	16 in.	16 in.
Overall Height:		
Overall Open Width:	25 in.	27 in.
Veight:	40 lbs.	44 lbs.
Veight Capacity:	450 lbs.	450 lbs.
ICPCS Code:	E1038	E1038

KM-TV10B Transport Chair







The Karman KM-TV10B travel chair is simply the heavy weight champion of all super lightweight folding travel chairs on the market today. Engineered to last, and affordable! It has a reinforced patented all 6061 T-6 aircraft-grade aluminum frame; the Karman travel chair offers maximum convenience, comfort, and security.





Color: Black

Champagne

- Padded flip back armrests
- Folds easily in half w/ quick release
- Fold up footrests
- Fits virtually in every car trunk
- Very light weight, easy to carry

Product Specifications	KM-TV-10B16	KM-TV-10B18
Seat Width:	16 in.	16 in.
Seat Depth:	18 in.	18 in.
Overall Open Width:		
Seat Height:	21 in.	21 in.
Back Height:	19 in.	19 in.
Overall Height:	36 in.	36 in.
Overall Width:	23.5 in.	25.5 in.
Weight:	18 lbs.	18 lbs.
Weight Capacity:	225 lbs.	225 lbs.
HCPCS Code:	E1038	E1038

MVP-502 Transport Chair



Product Features

The Karman Healthcare MVP-502 Ergonomic Transport Recliner offers "V" Shape seating that prevents sliding, sores, and relieves pressure from the body. Enjoy the reclining comfort with a backrest that effortlessly reclines from 93-160 degrees and ergonomic armrests that adjust when reclined.

Color: Black

- Folds ultra compact
- Backrest reclines 93-160 degrees
- Extending brakes
- Footrests stretch to reduce pressure
- Adjustable length anti-slip footplates

Product Specifications	MVP-502TP-16	MVP-502TP-18	MVP-502TP-20
Seat Width:	16 in.	18 in.	20 in.
Seat Depth:	18 in.	18 in.	18 in.
Armrest Height:			
Seat Height:	20 in.	20 in.	20 in.
Back Height:	19 in.	19 in.	19 in.
Overall Height:	52 in.	52 in.	52 in.
Overall Open Width:	23 in.	25 in.	27 in.
Folded Width:	14 in.	15 in.	16 in.
Weight:	34 lbs.	34 lbs.	34 lbs.
Weight Capacity:	250 lbs.	250 lbs.	250 lbs.
HCPCS Code:	K0003	K0003	K0003

VIP-515 Transport Chair





Product Features

The Karman Lightweight Tilt-in-Space VIP-515 is a revolutionary tilt chair that weighs just 34 pounds in transport.

This type of a chair redistributes pressure to the lower back as well as under the legs and it has a 0-35° weight-shifting center-of-gravity adjustment.

Color: Black

- Shifts pressure to reduce bedsores
- Flip back height-adjustable armrests
- Seat belt restraint included
- Center of gravity weight shifting
- Exclusive perfectly matched gas cylinders

Product Specifications	VIP-515-TP16	VIP-515-TP18
Seat Width:	16 in.	18 in.
Seat Depth:	17 in.	17 in.
Armrest Height:		
Seat Height:	20 in.	20 in.
Back Height:	18+ in.	18+ in.
Overall Length:	44 in.	44 in.
Folded Width:	15 in.	15 in.
Weight:	33 lbs.	33 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	N/A	N/A

KM-5000 Transport Chair




Product Features

Karma Ultra Lightweight Recliner Wheelchair KM-5000 is a high strength, lightweight wheelchair is in a league all its own. The reclining range is 90° to 163° and the compact folding design make it ideal for storage or transportation. Additionally, the lightweight, aircraft-grade aluminum frame makes the chair weigh only 36 lbs.

Color: Black

- Adjustable height pillow
- Adjustable toggle wheel locks
- Padded detachable leg rests
- Black powder coat finish
- High Strength Aluminum Frame

Product Specifications	KM-5000TP-16	KM-5000TP-18	KM-5000TP-20	KM-5000TP-2
Seat Width:	16 in.	18 in.	20 in.	22 in.
Seat Depth:	18 in.	18 in.	18 in.	18 in.
Seat Height:	20 in.	20 in.	20 in.	20 in.
Back Height:	33 in.	33 in.	33 in.	33 in.
Overall Height:	51 in.	51 in.	51 in.	51 in.
Overall Open Width:	25 in.	27 in.	29 in.	31 in.
Weight:	33 lbs.	33 lbs.	34 lbs.	34 lbs.
Weight Capacity:	250 lbs.	250 lbs.	300 lbs.	300 lbs.
HCPCS Code:	N/A	N/A	N/A	N/A



Stand Up Wheelchairs

XO-101 - 47

XO-202 - 48

XO-101 Stand Up Chair



Product Features

The XO-101 is the lightest power stand up wheelchair on the market, this wheelchair weighs 56 lbs. total w/ battery.

This wheelchair features padded flip back armrests, flip up footplates, and quick release detachable rear wheels, which comes in handy. Also available in FULL POWER STAND XO-202.

Color: Black

- Padded flip back armrests
- Flip up footplates
- Easy push button for stand up position
- Quick release rear wheels
- Stand up freeze prevention design

Seat Width:	16 in.	18 in.
Seat Depth:	18 in.	18 in.
Back Height:	16 - 19 in.	16 - 19 in.
Overall Height:		
Overall Open Width:	25 in.	25 in.
Weight:	56 lbs.	56 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	N/A	N/A

XO-202 Stand Up Chair



Product Features

The Karman Power Wheelchair with Power Stand-Up offers many health benefits previously unavailable to those in a wheelchair.

By putting the body in a standing position, the chair provides users with better blood circulation, kidney and bladder functioning, as well as muscle tone.

Color: Silver

- Adjustable footrests
- Adjustable backrest height
- Easy to use remote for drive & stand
- Flat free front and rear wheels
- See through sleek side panels

roduct Specifications	XO-202	XO-202N
Seat Width:	18 in.	16 in.
Seat Depth:	18 in.	18 in.
Back Height:	16 - 19 in.	16 - 19 in
Overall Height:		
Overall Open Width:	26.5 in.	25 in.
Weight:	110 lbs.	110 lbs.
Weight Capacity:	250 lbs.	250 lbs.
HCPCS Code:	N/A	N/A



Active Wheelchairs

AT-60 - **50**

S-ERGO ATX - 51

AT-60 Active Wheelchair



Product Features

The Karman AT-60 features a high-strength, rigid aircraft grade alloy frame, making it both lightweight and sturdy.

The backrest is foldable and angle adjustable, making the chair ideal for transport while allowing the user to find the perfect position to suit their needs.

Color: Black

- Aircraft grade alloy frame
- Retightening straps
- Folding backrest
- Quick release rear wheels
- Double adjustable handbrakes

Product Specifications	AT-60
Seat Width:	16 in.
Seat Depth:	17 in.
Back Height:	11 in.
Overall Height:	
Overall Open Width:	24 in.
Weight:	23 lbs.
Weight Capacity:	220 lbs.
HCPCS Code:	N/A

ATX Active Wheelchair



Product Features

The new ATX is ideal for all those who want to ride through life boldly. Thanks to outstanding agility and rigidity of the unique multi-geometrical extruded cross bars.

Special features of this "extreme design" include a compact dimensions with rigid structure.

Color: Green

- One piece foldable footplate
- Flip back armrest
- Seat height adjustable
- Tension adjustable back rest
- 15.4 lbs. starting frame weight

Product Specifications	S-ERGO-ATX14	S-ERGO-ATX16	S-ERGO-ATX18
Seat Width:	14 in.	16 in.	18 in.
Seat Depth:	15 in.	16 in.	18 in.
Seat Height:	14 in.	14 in.	14 in.
Back Height:			
Overall Height:	var.	var.	var.
Overall Open Width:	var.	var.	var.
Weight:	var.	15.4 lbs.	var.
Weight Capacity:	275 lbs.	275 lbs.	275 lbs.
HCPCS Code:	N/A	N/A	N/A



Rollators

R-3600 / R-4100 - 53

R-4100N / R-4200 - 54

R-4600 / R-4602-T - 55

R-4608 / R-4700 - 56

R-4800 / KW-100 - 57

R-3600 / R-4100 Rollator







Color: Burgundy Color: Blue

R-3600 Product Features

The Karman 3-Wheeler is a light 13-lbs. rollator/walker. Allows user to navigate obstacles and doorways in-passable with a standard 4-leg Walker. Three-Wheel rollators work great for people that just need something to keep their balance.

Includes both a basket and storage bag to keep everything you need.

R-4100 Product Features

The Karman R-4100 Series Rollator is a highly maneuverable fall prevention aid designed with a lower seat for individuals requiring a compact and lightweight personal mobility aid. Lighter than most two-wheel walkers at only 11 pounds, this durable and reliable rollator delivers the convenience you want at a price you can't pass up.

Veight: Veight Capacity:	13 lbs.	11 lbs.
Weight Capacity:		
	250 lbs.	250 lbs.
Handle Gap:	N/A	18 in.
Overall Height:	34.5 - 40.5 in.	28 - 34 in.
Brakes:	Loop Brake	Loop Brake
HCPCS Code:	E0143	E0143

R-4100N / R-4200 Rollator







Color: Burgundy Color: Blue

R-4100N Product Features

The Karman R-4100N Series Rollator is a highly maneuverable fall prevention aid designed with a lower seat for individuals requiring a compact and lightweight personal mobility aid. Lighter than most two-wheel walkers at only 11 pounds, this durable and reliable rollator delivers the convenience you want at a price you can't pass up.

R-4200 Product Features

The R-4200 rollator with push-down brakes from Karman healthcare is a lightweight rollator with an unique braking system. Apply light downward pressure to the rollator's handles to apply the brakes. This style of braking is perfect for those with hands that cannot squeeze traditional handle brakes.

Product Specifications	R-4100N	R-4200	
Weight:	11 lbs.	11 lbs.	
Weight Capacity:	250 lbs.	250 lbs.	
Handle Gap:	N/A	N/A	
Overall Height:	32 - 35 in.	34.5 - 40.5 in.	
Brakes:	Loop Brake	Loop Brake	
Backrest Height:	10 in.	9.5 in.	
Seat Width:	16 in.	18 in.	
HCPCS Code:	E0143	E0143	

R-4600 / R-4602 Rollator







Color: Burgundy

Color: Blue

R-4600 Product Features

The Karman 3-Wheeler is a light 13-lbs. rollator/walker. Allows user to navigate obstacles and doorways in-passable with a standard 4-leg Walker. Three-Wheel rollators work great for people that just need something to keep their balance.

Includes both a basket and storage bag to keep everything you need.

R-4602 Product Features

The Karman 4-wheel rollator and transport chair combines the functions of a rollator and a transport wheelchair into one lightweight device. The rollator turns into a transport chair by simply flipping down the footrests. This 4-wheel rollator/transport chair is very lightweight and folds easily for convenient storage and transport.

Product Specifications	R-4600	R-4602	
Weight:	11 lbs.	16 lbs.	
Weight Capacity:	250 lbs.	250 lbs.	
Handle Gap:	N/A	18 in.	
Overall Height:	34.5 - 40.5 in.	35 - 39 in.	
Brakes:	Loop Brake	Loop Brake	
HCPCS Code:	E0143	E0143	
Horos code.	E0143	E0143	

R-4608 / R-4700 Rollator







Color: Burgundy Color: Blue

R-4608 Product Features

This lightweight aluminum walker is great for travelling; it weighs only 15 lbs and can easily fold to fit in the trunk of your car.

The 8" and the loop brakes are easy to use while providing peace of mind. Functional loop brakes can also lock in place to make it easy to rest on the rollator.

R-4700 Product Features

R-4700w Extra Wide Steel Rollator is a heavy-duty rollator that features a wide 23" handle opening and a 400 lbs. weight capacity.

The steel frame is sturdy and strong and a padded flip-down seat allows the user to stop and rest anywhere they choose without a problem.

Product Specifications	R-4608	R-4700	
Weight:	15 lbs.	24 lbs.	
Weight Capacity:	250 lbs.	400 lbs.	
Handle Gap:	N/A	N/A	
Overall Height:	34 - 38 in.	33 - 37 in.	
Brakes:	Loop Brake	Loop Brake	
Backrest Height:	11 in.	10.5 in.	
Seat Width:	18 in.	23 in.	
HCPCS Code:	E0143	E0149	

R-4800 / KW-100 Walkers



Color: Burgundy Color: Blue Color: Burgundy Color: Black Color: White

R-4800 Product Features

This aluminum rollator by Karman is extra wide and features a 400 pound weight capacity for bariatric users. The rollator is super lightweight at just 20 pounds and easily folds for storage or transport.

KW-100 Product Features

The Knee Scooter offers increased comfort and ease of mobility to patients who cannot put weight on their foot or ankle, but want to remain active. With its four large wheels and padded knee platform, the Knee Scooter is perfect for indoor and outdoor use. The offset post on the aluminum knee platform allows for right-to-left adjustment in seconds - no tools required.

roduct Specifications	R-4800	KW-100	
Weight:	11 lbs.	22.5 lbs.	
Weight Capacity:	250 lbs.	250 lbs.	
Pad Size:	N/A	15 - 7 in.	
Overall Height:	34.5 - 40.5 in.	31 - 41 in.	
Overall Width:	N/A	16.5 in.	
Brakes:	Loop Brake	Loop Brake	
HCPCS Code:	E0143	E0118	

Daily Aids - Accesories



SC-200 Shower Chairs

- Easy assembly
- Blow molded
- 3 adjustable positions
- Lightweight aluminum frame
- Angled legs with suction style tips



SC-100 Shower Chairs

- Easy assembly
- Drainage holes
- 3 adjustable positions
- Lightweight aluminum frame
- Angled legs with suction style tips



PE2 Exercise Equipment

- Stimulates circulation
- tension adjustable
- Box assembled
- Compact and durable



OXYGEN REGULATORS

- Reduces pressure
- Controls flow of oxygen
- Aluminum allow body
- Easy to read flow rates
- All brass core



OT-10 Overbed Table

- Chrome plated steel
- Weight capacity 40 lbs.
- Swivel casters
- Table top can be raised
- Overall width at the base



NEBULIZERS

- Durable piston
- Fast, Effective treatments
- Medical capacity of 5 ml.
- Lightweight 4.3 lbs.



TC-101 Tire Changer

- Flexible design
- Fits 4 26 in.
- Easy to use
- Durable



UNIVERSAL Folding Headrest

- Fits 16 18 in. seats
- Fits 20 22 in. seats
- Folds with wheelchair
- Great durability



Luxury Canopy

- Easy assembly
- Folds with wheelchair
- Adjustable positions
- Can rotate back to front
- Fits 16 20 in, back width



www.karmanhealthcare.com 19255 San Jose Ave. City of Industry, CA 91748 Phone 626-581-2235 Fax 626-581-2335 © 2013 Karman Healthcare

Karma



PART#

KARMAN Healthcare Inc.

19255 San Jose Ave.

DESCRIPTION

City of Industry, CA 91748 • Tel: 626-581-2235 • Fax: 626-581-2335

SUGGESTED

Wheelchair Parts and Accessories Price List

1AKI#	<u>BESOKII ITON</u>	RETAIL PRICE
ANTI-TIPPERS AT105	Anti-tipper for S-Ergo 105 (pr)	Pair \$43
AT115	Anti-tipper for S-Ergo 115 (pr)	\$43
AT305	Anti-tipper for S-Ergo 305 (pr)	\$50
AT800-DY	Anti-tipper, 7/8" Tube for KN-800B, LT-800B, LT-700B (pr)	\$32
AT44	Anti-tipper, Pin On, For LT-K4 (pr)	\$32
AT802	Anti-tipper, Pin On, For KM802 (pr)	\$43
AT802DY	Anti-tipper, Pin On, For 802-DY (pr)	\$32
AT8020	Anti-tipper, Pin On, For KM8000, KM5000 (pr)	\$43
AT9020L	Anti-tipper, Pin On, For KM-9020L (pr)	\$32
AT950	Anti-tipper, Pin On, For LT-950 (pr)	\$32
ARMRESTS ASSE	MBLIES	Each
AR22L	Desk Length, Side Panel, Left	\$50
AR22R	Desk Length, Side Panell, Right	\$50
AR23L	Full Length, Side Panel, Left	\$50
AR23R	Full Length, Side Panel, Right	\$50
AR87L	Desk Swing Back, Aluminum Panel, KM8000, Left	\$68
AR87R	Desk Swing Back, Aluminum Panel, KM8000, Right	\$68
AR88L	Full, Swing Back, Aluminum Panel, KM8000, Left	\$68
AR88R	Full, Swing Back, Aluminum Panel, KM8000, Right	\$68
ARMREST PADS		Each
AP21	Grey Plastic, Universal Full Arm	\$9
AP22B	Upholstered, Universal Full Arm Pad, Black	\$9
AP22U	Upholstered, Universal Full Arm Pad, Dark Blue	\$9
AP23B	Upholstered, Universal Desk Arm Pad, Black	\$9
AP23U	Upholstered, Universal Desk Arm Pad, Dark Blue	\$9
AP87	Upholstered, Desk Arm Pad, Black for KM-8000	\$16
AP88	Upholstered, Full Arm Pad, Black for KM-8000	\$16
AP-105	Molded Ureform, Wide Arm Pad, Black for S-Ergo 105/115	\$29
AP-305	Molded Ureform, Wide Arm Pad, Black for S-Ergo 305	\$29
BASKET		Each
BASKET-14	Basket 14.5"x10"x 6" Small for R-4100N/ R-4200/ R-4600old style	\$16
BASKET-15	Basket 15.5"x10"x 6" Small for R-4100/ R-4200 /R-4600	\$16
BASKET-3600	Basket for R-3600	\$16
BAG		Each
BAG-TV10A	Bag for TV-10A	\$22
BRAKE LOCK ASS		Each
BR700L	Push Brakes, Detachable Arm, Tube Mount, Left	\$16
BR700R	Push Brakes, Detachable Arm, Tube Mount, Right	\$16
BR800L	Push Brakes, Fixed Arm, Tube Mount, Left	\$16
BR800R	Push Brakes, Fixed Arm, Tube Mount, Right	\$16
BR-950L	Push Brakes for LT-950 Left Side	\$22
BR-950R	Push Brakes for LT-950 Right Side	\$22
BR1000L	Push Brakes, Transporter, Tube Mount, Left	\$16
BR1000R	Push Brakes, Transporter, Tube Mount, Right	\$16
BR8000L	Push Brakes, KM-8000 Lt Wt, Tube Mount, Left	\$34

Page 1

BR8000R BR5000L BR5000R BR4600	Push Brakes, KM-8000 Lt Wt, Tube Mount, Right Push Brakes, KM-5000 Lt Wt, Tube Mount, Left Push Brakes, KM-5000 Lt Wt, Tube Mount, Right Brakes Assembly w/Handle for R-4600	\$34 \$34 \$34 \$25	
BRAKE EXTENSIO BE6 BE8 BE11 BE88	N Brake Extension, 5/8"space x 6" length Brake Extension, 5/8" space x 8" length Brake Extension, 5/8" space x 11" length Brake Extension, for KM8000	Each \$9 \$9 \$9 \$13	
BRAKE TIPS BT22 BT88	Brake Tips , Black, Flat Brake Tips , Black for KM8000	Each \$4 \$4	
BUMPERS B7 B8	Bumper Black, 7/8" Bumper Black, 1"	Each \$4 \$4	
CASTERS w/Bearin C42B C46B C81B C81G C51LT C61LT C71LT C81LT CW3 CW5	Front Caster for R-4200 Front Caster for R-4600 8" x 1" Front Caster 8" x 1" Front Caster, Grey Tire 5" x 1" Light Weight Caster, w/ Polyurethane 6" x 1" Light Weight Caster, w/ Polyurethane 7" x 1" Light Weight Caster, w/ Polyurethane 8" x 1" Light Weight Caster, w/ Polyurethane 8" x 1" Light Weight Caster, w/ Polyurethane 3" Wheel w/Tube for W-01-W3 (pr) 5" Wheel w/Tube for W-01-W5 (pr)	Each \$14 \$14 \$20 \$20 \$25 \$25 \$25 \$25 \$25 \$25 \$25	
CASTER FORKS CF22 CF88 CF99	Caster Fork Standard Steel Caster Fork, Light Weight for KM8000 Caster Fork, Heavy Duty	Each \$22 \$27 \$27	
CROSS BRACES CB2 CB2N CB88 CB88N	Cross Braces 18" Width Cross Braces 16" Width Cross Braces 18" Width for Lt Weight Cross Braces 16" Width for Lt Weight	Pair \$76 \$76 \$86 \$86	
CUSHION CU16 CU18 CU-XO202-SEAT CU-XO202N-SEAT CU-XO202-BACK CU-XO202-BACK	High Density 2" Seat Cushion, 16"x16" High Density 2" Seat Cushion, 18"x16" Memory Seat 2" Cushion 18"x17" Memory Seat 2" Cushion 16"x17" Memory Cushion Back for 18" Seat Width Memory Cushion Back for 16" Seat Width	Each \$68 \$68 \$158 \$158 \$86 \$86	
ELEVATING LEGR EL18BB-INV-DY EL16BB-INV-DY EL18CB-INV-DY EL16CB-INV-DY EL18CB-EJ-DY EL8000 EL-S305R-18 EL-S305S-18		Pair \$97 \$97 \$97 \$97 \$97 \$228 \$228 \$228 \$228	Page 2

EL-S305S-16 EL-VIP515-18 EL-VIP515-16	Elevating Legrests, Padded, for S-ERGO-305 16" Silver (pr) Elevating Legrests, Padded, for VIP-515 18" Black Diamond (pr) Elevating Legrests, Padded, for VIP-515 16" Black Diamond (pr)	\$228 \$228 \$228
REFLECTION REA RF-115	R WHEEL COVERS Reflection Rear Wheel Cover for	Pair \$65
SWING-AWAY FOO FR14 FR16 FR18 FR18B-INV-DY FR8000 FR-S115R-18 FR-S115S-16 FR-S115S-16 FR-S305R-16 FR-S305R-16 FR-S305S-18 FR-S305S-18	Swing-away Footrests, 14" Chairs (pr) Swing-away Footrests, 16" Chairs (pr) Swing-away Footrests, 18" to 24" Chairs, Also for Transporter (pr) Swing-away Footrests, Black Frame (pr) Swing-away Footrests, for KM8000 Lt Wt, (pr) Swing-away Footrests, for S-ERGO115 18" Red (pr) Swing-away Footrests, for S-ERGO115 16" Red (pr) Swing-away Footrests, for S-ERGO115 18" Silver (pr) Swing-away Footrests, for S-ERGO115 16" Silver (pr) Swing-away Footrests, for S-ERGO305 18" Red (pr) Swing-away Footrests, for S-ERGO305 16" Red (pr) Swing-away Footrests, for S-ERGO305 18" Silver (pr) Swing-away Footrests, for S-ERGO305 16" Silver (pr)	Pair \$62 \$62 \$62 \$126 \$126 \$126 \$126 \$126 \$1
FR-VIP515-18 FR-VIP515-16	Swing-away Footrests, for VIP-515 18" Black Diamond (pr) Swing-away Footrests, for VIP-515 16" Black Diamond (pr)	\$126 \$126
HAND GRIPS HG6 HG7 HG8 HG46	Hand Grips 3/4" Tube, Black, for KM8000 Lt Wt Hand Grips 7/8" Tube, Black Hand Grips 1" Tube, Black, for KM8000 Lt Wt Hand Forms for R-4200/R-4600	Each \$8 \$6 \$8 \$4
HANDLES HD-LT2019 HD-R4600 HD-XO202	Handle -LT2019 Handle with Loop Brake Cable- R-4600 Push Handle for XO101/XO202	Each \$25 \$25 \$68
HAND RIMS HR22 HR8000	Black, Composite Hand Rim for KN800 Black, Composite Hand Rim for KM8020	Each \$32 \$32
HEADREST HDRST-FLD-115	Foldable Headrest Height Adjustable for 7/8" Frame	Each \$196
HEEL LOOPS HL18	Heel Loops	Each \$9
I.V. POD ATTACHI IV2	I.V. Pod Attachment	Each \$45
LEGREST PADS LP5B LP5U LP8	Upholstered, Padded Calf Support, Black Upholstered, Padded Calf Support, Dark Blue Upholstered, Padded Calf Support, Black KM8000	Each \$18 \$18 \$18
LEG STRAPS LS16 LS18	Leg Strap, Black Leg Strap, Black	Each \$11 \$11

OXYGEN TANK H	OLDER	Each
OX2	Oxygen Tank Holder for Standard Wheelchair	\$50
OX88	Oxygen Tank Holder for KM-8000	\$54
REAR WHEEL AS	SEMBLIES	Each
RW12	Rear Mag Wheel Black, 12 1/2" x 2 1/4", for T-1000	\$45
RW12PU	Rear Mag Wheel Black, 12 1/2" x 2 1/4" PU, for T-1000PU	\$50
RW241	Rear Mag Wheel Black, 24" x 1" Rubber Tire, for KN Steel Series	\$86
RW44	Rear Mag Wheel Black, 24" x 1", for LT-K4	\$90
RW88	Rear Mag Wheel Black, 24" x 1" PU Tire, for KM8000	\$122
RW-S105	Rear Spoke Wheel, 24" x 1" PU Tire, Fixed for S-Ergo105/115	\$122
RW-S305Q	Rear Spoke Wheel, 24" x 1" PU Tire, Quick Release for S-Ergo305	\$122
RECLINING MECH		Each
RM22L	Left Mechanism for KN-880	\$68
RM22R	Right Mechanism for KN-880	\$68
RM88L	Left Mechanism for KM-5000	\$86
RM88R	Right Mechanism for KM-5000	\$86
RM-VIP515-L	Left Mechanism for VIP-515	\$104
RM-VIP515-L	Left Mechanism for VIP-515	\$104
SEAT GUIDES		Each
SG800	Seat Guides Plastic, Fixed Arms	\$4
SG700	Seat Guides Plastic, Detachable Arms	\$4
SEAT BELTS		Each
SB22	Seat BeltTwo-Piece Seat Belt. 45" Long, 2" Webbing.	\$14
SB88	One-Piece Seat Belt. 60" Long, 2" Webbing.	\$14
SB99	Seat Belt - Push Button, Auto Style. 60" Long, 2" Webbing.	\$20
SIDE DANEI S		Each
SIDE PANELS	Sido Panol for S. Ergo 205 Loft Sido	Each
SP-S305-L	Side Panel for S-Ergo 305 Left Side	\$29
SP-S305-L SP-S305-R	Side Panel for S-Ergo 305 Right Side	
SP-S305-L SP-S305-R TIRES AND TUBE	Side Panel for S-Ergo 305 Right Side	\$29 \$29 Each
SP-S305-L SP-S305-R TIRES AND TUBE T8141	Side Panel for S-Ergo 305 Right Side S Pneumatic Tire 8" x 1 1/4"	\$29 \$29 Each \$10
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241	Side Panel for S-Ergo 305 Right Side S Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4"	\$29 \$29 Each \$10 \$13
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8"	\$29 \$29 Each \$10 \$13 \$13
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1"	\$29 \$29 Each \$10 \$13 \$13 \$14
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8"	\$29 \$29 Each \$10 \$13 \$13 \$14 \$13
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1"	\$29 \$29 Each \$10 \$13 \$13 \$14 \$13 \$29
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1"	\$29 \$29 Each \$10 \$13 \$13 \$14 \$13 \$29 \$29
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Inner Tube 8" x 1 1/4"	\$29 \$29 Each \$10 \$13 \$13 \$14 \$13 \$29 \$29 \$9
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4"	\$29 \$29 Each \$10 \$13 \$14 \$13 \$29 \$29 \$9
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU8141 TU12241 TU22183	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8"	\$29 \$29 Each \$10 \$13 \$14 \$13 \$29 \$29 \$9 \$9
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241 TU22183 TU24183	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8" Inner Tube 24" x 1 3/8"	\$29 \$29 Each \$10 \$13 \$14 \$13 \$29 \$29 \$9
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU8141 TU12241 TU22183	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8"	\$29 \$29 Each \$10 \$13 \$14 \$13 \$29 \$29 \$9 \$9
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241 TU22183 TU24183	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8" Inner Tube 24" x 1 3/8"	\$29 \$29 Each \$10 \$13 \$13 \$14 \$13 \$29 \$29 \$9 \$9 \$9
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU8141 TU12241 TU22183 TU24183 TP241	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8" Inner Tube 24" x 1 3/8"	\$29 \$29 Each \$10 \$13 \$14 \$13 \$29 \$29 \$9 \$9 \$9
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241 TU22183 TU22183 TU24183 TP241 UPHOLSTERY	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8" Inner Tube 24" x 1 3/8" P/U Polyurethane Tire 24" x 1"	\$29 \$29 Each \$10 \$13 \$13 \$14 \$13 \$29 \$9 \$9 \$9 \$9
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241 TU22183 TU22183 TU24183 TP241 UPHOLSTERY US1416	Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8" Inner Tube 24" x 1 3/8" P/U Polyurethane Tire 24" x 1" Seat 14"W x 16"L Seat 16"W x 16"L	\$29 \$29 Each \$10 \$13 \$14 \$13 \$29 \$29 \$9 \$9 \$9 \$9 \$9 \$22 Each \$27
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241 TU22183 TU22183 TU24183 TP241 UPHOLSTERY US1416 US1616	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8" Inner Tube 24" x 1 3/8" P/U Polyurethane Tire 24" x 1" Seat 14"W x 16"L	\$29 \$29 Each \$10 \$13 \$14 \$13 \$29 \$29 \$9 \$9 \$9 \$9 \$9 \$22 Each \$27 \$27
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241 TU22183 TU24183 TP241 UPHOLSTERY US1416 US1616LT	S Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8" Inner Tube 24" x 1 3/8" P/U Polyurethane Tire 24" x 1" Seat 14"W x 16"L Seat 16"W x 16"L Light Weight for KM-8000 Seat 18"W x 16"L	\$29 \$29 Each \$10 \$13 \$13 \$14 \$13 \$29 \$29 \$9 \$9 \$9 \$9 \$9 \$22 Each \$27 \$41
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241 TU22183 TU24183 TP241 UPHOLSTERY US1416 US1616 US1616LT US1816	Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 1/4" Pneumatic Tire 24" x 1 1/4" Airless Insert 22" x 1 1/4" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8" Inner Tube 24" x 1 3/8" P/U Polyurethane Tire 24" x 1" Seat 14"W x 16"L Seat 16"W x 16"L Seat 16"W x 16"L Light Weight for KM-8000	\$29 \$29 Each \$10 \$13 \$13 \$14 \$13 \$29 \$9 \$9 \$9 \$9 \$9 \$22 Each \$27 \$27 \$41 \$27 \$41
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241 TU22183 TU24183 TP241 UPHOLSTERY US1416 US1616 US1616LT US1816 US1816LT US1818LT	S Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8" Inner Tube 22" x 1 3/8" Inner Tube 24" x 1 3/8" P/U Polyurethane Tire 24" x 1" Seat 14"W x 16"L Seat 16"W x 16"L Seat 18"W x 16"L Light Weight for KM-8000	\$29 \$29 Each \$10 \$13 \$14 \$13 \$29 \$29 \$9 \$9 \$9 \$9 \$22 Each \$27 \$27 \$41 \$27 \$41 \$41
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241 TU22183 TU22183 TU24183 TP241 UPHOLSTERY US1416 US1616 US1616LT US1816 US1816LT US1818LT US2016	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 24" x 1 3/8" P/U Polyurethane Tire 24" x 1" Seat 14"W x 16"L Seat 16"W x 16"L Seat 16"W x 16"L Seat 18"W x 16"L	\$29 \$29 Each \$10 \$13 \$14 \$13 \$29 \$29 \$9 \$9 \$9 \$9 \$22 Each \$27 \$27 \$41 \$27 \$41 \$41 \$40
SP-S305-L SP-S305-R TIRES AND TUBE T8141 T12241 T22183 T241 T24183 T1221 T1241 TU8141 TU12241 TU22183 TU24183 TP241 UPHOLSTERY US1416 US1616 US1616LT US1816 US1816LT US1818LT	Side Panel for S-Ergo 305 Right Side Pneumatic Tire 8" x 1 1/4" Pneumatic Tire 12 1/2" x 2 1/4" Pneumatic Tire 22" x 1 3/8" Rubber Tire 24" x 1" Pneumatic Tire 24" x 1 3/8" Airless Insert 22" x 1" Airless Insert 22" x 1" Airless Insert 24" x 1" Inner Tube 8" x 1 1/4" Inner Tube 12 1/2" x 2 1/4" Inner Tube 22" x 1 3/8" Inner Tube 24" x 1 3/8" P/U Polyurethane Tire 24" x 1" Seat 14"W x 16"L Seat 16"W x 16"L Seat 16"W x 16"L Seat 18"W x 18"L for KM-5000	\$29 \$29 Each \$10 \$13 \$14 \$13 \$29 \$29 \$9 \$9 \$9 \$9 \$22 Each \$27 \$27 \$41 \$27 \$41 \$41

Page 4

US2218 US2218LT US2418 US2618	Seat 22"W x 18"L Seat 22"W x18"L Light Weight for KM-8000 Seat 24"W x 18"L Seat 26"W x 18"L	\$50 \$65 \$54 \$63
UB1418 UB1618 UB1617LT UB1818 UB1817LT UB2018 UB2018LT UB2218 UB2218LT UB2418 UB2618	Back 14"W x 18"L Back 16"W x 18"L Back 16"W x 17"L Light Weight for KM-8000 Back 18"W x 18"L Back 18"W x 17"L Light Weight for KM-8000 Back 20"W x 18"L Back 20"W x17"L Light Weight for KM-8000 Back 22"W x 18"L Back 22"W x 18"L Back 22"W x 18"L Back 24"W x 18"L Back 26"W x 18"L	\$27 \$27 \$41 \$27 \$41 \$45 \$59 \$50 \$65 \$54
UPPER UPHOLST	ERY REMOVABLE	Each
US-S115 US-S115N US-S115-20 US-S305 US-S305N	Removable Upper Seat 18"x17" for S-Ergo 105, 106 & 115 Removable Upper Seat 16"x17" for S-Ergo 105, 106 & 115 Removable Upper Seat 20"x17" for S-Ergo 105 & 115 Removable Upper Seat 18"x17" for S-Ergo 305 Removable Upper Seat 16"x17" for S-Ergo 305	\$50 \$50 \$50 \$50 \$50 \$50
UB-S115 UB-S115N UB-S115-20 UB-S305 UB-S305N	Removable Upper Back 18"x17" for S-Ergo 105 & 115 Removable Upper Back 16"x17" for S-Ergo 105 & 115 Removable Upper Back 20"x17" for S-Ergo 115 Removable Upper Back 18" for S-Ergo 305 Removable Upper Back 16" for S-Ergo 305	\$50 \$50 \$50 \$50 \$50