

karma



Lightweight

Wheelchairs

LT-700 Lightweight Class

STANDARD LIGHTWEIGHT WITH DETACHABLE ARM

SPECIFICATIONS	LT-700NT	LT-700T
Seat Width	16 in.	18 in.
Seat Depth	16 in.	16 in.
Armrest Height	10 in.	10 in.
Seat Height	17/19 in.	17/19 in.
Back Height	17 in.	17 in.
Overall Width	36 in.	36 in.
Overall Open Width	23 in.	25 in.
Folded Width	11 in.	11 in.
Overall Length	31 in.	31 in.
Weight w/o Riggings	36 lbs.	36 lbs.
Weight Capacity	250 lbs.	250 lbs.
HCPSC Code	K0003	K0003

(SUBJECT TO CHANGE WITHOUT NOTICE)



KEEPING **YOU** IN MIND...

- Seat height adjustable / dual axle
- Detachable / Padded desk length armrests with allow side guard
- Swing away/ Length adjustable footrests with heel loop footplates
- Push to lock wheel brakes
- 24 x 1" Mag wheels with polyurethane tires
- Extra-durable forks
- 8 x 1 flat free front casters
- High-grade / Padded nylon flame retardant upholstery
- Seat width: 16" or 18"
- Frame color: Silver vein

LT-800 Lightweight Class

STANDARD LIGHTWEIGHT WITH DETACHABLE ARM

SPECIFICATIONS	LT-800NT	LT-800T
Seat Width	16 in.	18 in.
Seat Depth	16 in.	16 in.
Armrest Height	11 in.	11 in.
Seat Height	19 in.	19 in.
Back Height	16 in.	16 in.
Overall Width	36 in.	36 in.
Overall Open Width	22 in.	24 in.
Folded Width	11 in.	11 in.
Overall Length	31 in.	31 in.
Weight w/o Riggings	34 lbs	34 lbs.
Weight Capacity	250 lbs.	250 lbs.
HCPSC Code	K0003	K0003

(SUBJECT TO CHANGE WITHOUT NOTICE)



MADE FOR DURABILITY...

- Fixed padded full armrests with alloy side guard
- Swing-away detachable footrests with heel loop
- Push to lock wheel brakes
- 24" Flat-free polyurethane rear tires
- Extra-durable forks
- 8" Flat free polyurethane casters
- High-grade / Padded nylon flame retardant upholstery
- Seat width: 16" or 18"
- Frame color: Silver vein

Lightweight
Wheelchairs

LT-980 Lightweight Class

STANDARD LIGHTWEIGHT ALUMINUM HIGH STRENGTH

SPECIFICATIONS	LT-980	LT-980BD
Seat Width	18 in.	18 in.
Seat Depth	16 in.	16 in.
Armrest Height	9 in.	9 in.
Seat Height	19 in.	19 in.
Back Height	19 in.	19 in.
Overall Height	41 in.	41 in.
Overall Open Width	26 in.	26 in.
Folded Width	13 in.	13 in.
Overall Length	28 in.	28 in.
Weight w/o Riggings	24 lbs	24 lbs.
Weight Capacity	250 lbs.	250 lbs.
HCPSC Code	K0004	K0004

(SUBJECT TO CHANGE WITHOUT NOTICE)



BE IN **YOUR** BEST FORM...

Fixed armrests with molded arm pads

Folding backrest for travel or storage

Durable molded arm pads

24" Flat-free polyurethane rear tires

Breathable black nylon upholstery and black frame

6" x 1" Flat-free front polyurethane high profile wheels

Easy-grip contoured push-rims

Ultra lightweight 24 lbs w/o footrest

Seat width: 16" or 18"

Frame color: Silver vein

Lightweight Wheelchairs

karma



Karman Healthcare

19255 San Jose Ave.

City of Industry, CA

Telephone Number - 626-581-2235

www.karmanhealthcare.com