

Aluminum/Steel Foldable Rollators

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Seat-Hand	Seat Depth	Weight
R-3600	N/A	28"	N/A	26 1/2"	N/A	N/A	N/A	34 1/2"~40"	N/A	13 lbs
R-4100	18"	23"	N/A	21"	N/A	19"	13"	28~34"	13 1/2"	11 lbs
R-4100-N	16"	22 1/2"	N/A	18 3/4"	N/A	19"	13"	30"~34 1/2"	12"	11 lbs
R-4200	18"	23"	N/A	23"	N/A	22"	13 1/2"	33 1/2"~37"	14"	11 lbs
R-4600	18"	23 1/2"	N/A	23"	N/A	22"	9 1/2"	34"~39"	12"	11 lbs
R-4602-T	18"	24"	N/A	24 1/2"	24 1/2"	22"	8~13"	34 3/4"~39"	12"	16 lbs
R-4608	18"	24 1/2"	N/A	27 1/4"	N/A	22"	14"	34"~39 1/2"	14"	15 lbs
R-4700W	23"	28 1/2"	N/A	24 1/2"	N/A	22"	10 1/2"	33~37"	13 1/2"	24 lbs
R-4800W	23"	28 1/2"	N/A	24 1/2"	N/A	22"	10 1/2"	33~37"	13 1/2"	20 lbs
R-5200	18"	28 1/2"	N/A	23"	N/A	21"	11"	33~37"	13"	11 lbs
R-5600	18"	28 1/2"	N/A	23"	N/A	22"	11"	33~37"	13"	11 lbs

W Extra Wide Seat Width

E Elevating Legrest

N Narrow Seat Width

Transport Chairs

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Seat-Hand	Seat Depth	Weight
T-2017	17"	22 1/2"	11"	23"	34"	19"	19"	40"	16"	23 lbs.
T-2019	19"	23 1/2"	11"	23"	34"	19"	19"	40"	16"	23 lbs.
T-2700N	17"	21"	11"	23 1/2"	35"	20"	17"	37"	16"	29 lbs.
T-2700	19"	23"	11"	23 1/2"	35"	20"	17"	37"	16"	29 lbs.
LT-2017	17"	19"	10"	23"	34"	20"	19"	40"	16"	19 lbs.
LT-2019	19"	21"	10"	23"	34"	20"	19"	40"	16"	19 lbs.
LT-1000HB	19"	21 1/2"	10 1/2"	27 1/2"	37 1/2"	19"	19"	40"	16"	23 lbs.

HB Hand Brakes

N Narrow

N Narrow Seat Width

*Specifications are subject to change without notice.

Travel Chair

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Seat-Hand	Seat Depth	Weight
KM-TV10B18	18"	26"	26"	26"	37"	20"	20"	38"	17"	18 lbs
KM-TV10B16	16"	24"	26"	26"	37"	20"	20"	38"	17"	18 lbs

TV Travel Chair

Standard Wheelchair

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Seat-Hand	Seat Depth	Weight
KN-700T	18"	24"	11"	31"	42"	17"/19"	16"	36"	16"	39 lbs
KN-700NT	16"	26"	11"	31"	42"	17"/19"	16"	36"	16"	39 lbs
KN-800T	18"	25"	11"	31"	42"	19"	16"	36"	16"	37 lbs
KN-800NT	16"	23"	11"	31"	42"	19"	16"	36"	16"	37 lbs
LT-800T	18"	24"	11"	31"	42"	19"	16"	36"	16"	35 lbs.
LT-800NT	16"	22"	11"	31"	42"	19"	16"	36"	16"	35 lbs.
LT-700T	18"	25"	11"	31"	42"	17"/19"	16"	36"	16"	35 lbs.
LT-700NT	16"	23"	11"	31"	42"	17"/19"	16"	36"	16"	35 lbs.

N Narrow Seat Width

J Junior

*Specifications are subject to change without notice.

Bariatric Wheelchairs

MODEL No.										
KN-920W	20"	28"	13"	32 1/2"	44 1/2"	21"	18"	36"	16"/18"	41 lbs.
KN-920WE	20"	28"	13"	32 1/2"	44 1/2"	21"	18"	36"	16"/18"	41 lbs.
KN-922W	22"	30"	13 1/2"	32 1/2"	44 1/2"	22"	18"	36"	16"/18"	49 lbs.
KN-922WE	22"	30"	13 1/2"	32 1/2"	44 1/2"	22"	18"	36"	16"/18"	49 lbs.
KN-924W	24"	32"	14"	32 1/2"	44 1/2"	22"	18"	36"	18"	53 lbs.
KN-924WE	24"	32"	14"	32 1/2"	44 1/2"	22"	18"	36"	18"	53 lbs.
KN-926W	26"	34"	14 1/2"	32 1/2"	44 1/2"	22"	18"	36"	18"	58 lbs.
KN-926WE	26"	34"	14 1/2"	32 1/2"	44 1/2"	22"	18"	36"	18"	58 lbs.
KM-8520F-20W	20"	29"	16"	32"	42"	20"	15"	36"	18"	36 lbs.
KM-8520F-22W	22"	31"	16"	32"	42"	20"	15"	36"	18"	36 lbs.

W Extra Wide Seat Width

E Elevating Legrest

*Specifications are subject to change without notice.

Recliners

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Floor-Hand	Seat Depth	Weight
KN-880E	18"	26 1/2"	12 1/2"	35 1/2"	47 1/2"	20"	18"	35 1/2"	18"	50 lbs.
KN-880NE	16"	24 1/2"	11 1/2"	35 1/2"	47 1/2"	20"	18"	35 1/2"	18"	50 lbs.
KM-5000F16	16"	25"	15"	34 1/2"	47 1/2"	19 1/2"	19"	37"	18"	36 lbs.
KM-5000F18	18"	27"	15"	34 1/2"	47 1/2"	19 1/2"	19"	37"	18"	36 lbs.
KM-5000TP-16	16"	25"	13"	29 1/2"	42 1/2"	19 1/2"	19"	37"	18"	33 lbs.
KM-5000TP-18	18"	27"	13"	29 1/2"	42 1/2"	19 1/2"	19"	37"	18"	33 lbs.
KM-5000TP-20W	20"	28"	13 1/2"	30 1/2"	42 1/2"	19 1/2"	19"	37"	18"	35 lbs.
KM-5000TP-22W	22"	30"	13 1/2"	30 1/2"	42 1/2"	19 1/2"	19"	37"	18"	35 lbs.
KM-5000F-20W	20"	29"	16"	47 1/2"	42 1/2"	19 1/2"	19"	37"	18"	37 lbs.
KM-5000F-22W	22"	31"	16"	47 1/2"	42 1/2"	19 1/2"	19"	37"	18"	37 lbs.
MVP-502TP-16	16"	23"	14 1/2"	34 1/2"	45"	21 1/2"	21"	37"	19"	34 lbs.
MVP-502TP-18	18"	25"	14 1/2"	34 1/2"	45"	21 1/2"	21"	37"	19"	34 lbs.
MVP502TP-20	20"	27"	14 1/2"	34 1/2"	45"	21 1/2"	21"	37"	19"	34 lbs.
MVP502-16	16"	24"	14 1/2"	34 1/2"	45"	21 1/2"	21"	37"	19"	36 lbs.
MVP502-18	18"	26"	14 1/2"	34 1/2"	45"	21 1/2"	21"	37"	19"	36 lbs.
MVP502-20	20"	28"	14 1/2"	34 1/2"	45"	21 1/2"	21"	37"	19"	36 lbs.

N Narrow Seat Width

W Wide Seat Width

Ergonomic Wheelchair w/ Anti-Biotic Upholstery

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Floor-Hand	Seatdepth	Weight
S-105F16	16"	22.5"	12"	28"	39"	19"	15"	36"	17"	27 lbs
S-105F18	18"	25"	12"	28"	39"	19"	17"	36"	17"	27 lbs
S-106F16	16"	23"	12"	28"	39"	19"	17"	36"	17"	27 lbs.
S-106F18	18"	25"	12"	28"	39"	19"	17"	36"	17"	27 lbs
S-115F16	16"	22.5"	12"	28"	39"	19"	17"	36"	17"	25 lbs
S-115F18	18"	25"	12"	28"	39"	19"	17"	36"	17"	25 lbs
S-115F20	20"	27"	12"	28"	39"	19"	17"	36"	17"	26 lbs.
S-115Q16	16"	22.5"	12"	28"	39"	19"	17"	36"	17"	25 lbs
S-115Q18	18"	25"	12"	28"	39"	19"	17"	36"	17"	25 lbs
S-115Q20	18"	27"	12"	28"	39"	19"	17"	36"	17"	26 lbs
S-115MAG18	18"	25"	12"	28"	39"	19"	16"	36"	17"	26 lbs
S-115MAG16	16"	23"	12"	28"	39"	19"	16"	36"	17"	26 lbs.
S-115TP-16	16"	22 1/2"	12"	25"	35 3/4"	19"	16"	36"	17"	22 lbs
S-115TP-18	18"	23 1/2"	12"	25"	35 3/4"	19"	16"	36"	17"	22 lbs
S-115TP-20	20"	25"	12"	25"	35 3/4"	19"	16"	36"	17"	22 lbs.

Ergonomic Wheelchair w/ Anti-Biotic Upholstery (CONT)

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Floor-Handl	Seatdepth	Weight
S115TPWB16	16"	22 1/2"	12"	25"	35 3/4"	19"	16"	36"	17"	23 lbs.
S115TPWB18	18"	23 1/2"	12"	25"	35 3/4"	19"	16"	36"	17"	23 lbs.
S-2501TP-16	16"	21 3/4"	12"	24"	35 3/4"	18 1/2"	17"	36"	17"	18 lbs
S-2501TP-18	18"	23 3/4"	12"	24"	35 3/4"	18 1/2"	17"	36"	17"	18 lbs
S-305Q16	16"	22.5"	12"	30"	43"	19"	17"	38"	17"	29 lbs
S-305Q18	18"	25"	12"	30"	43"	19"	17"	38"	17"	29 lbs

Q Quick Release Wheels WB Wire Brakes F Fixed Armrest

*Specifications are subject to change without notice.

Ultra Lightweight Aluminum Wheelchair

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Floor-Hand	Seat Depth	Weight
LT-980	18"	26"	13 1/2"	29 1/2"	37"	19"	19"	39"	15"	25 lbs
KM-9020L	18"	25 1/4"	13"	26"	35 1/2"	19"	19"	37"	15"	25 lbs

*Specifications are subject to change without notice.

Tilt-in-Space Wheelchair & Transport

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Floor-Hand	Seat Depth	Weight
VIP515TP-16	16"	24"	15"	29"	39"	19"	18"	38"	16 1/2"	34 lbs
VIP515TP-18	18"	27"	15"	29"	39"	19"	18"	38"	16 1/2"	34 lbs
VIP515-18 MS	18"	27 1/2"	16"	31 1/2"	41"	19"	18"	39"	16 1/2"	38 lbs
VIP515-16 MS	16"	26"	16"	31 1/2"	41"	19"	18"	39"	16 1/2"	38 lbs.

TP Transport MS 20" Rear Wheels

High Strength Aluminum Wheelchairs

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Floor-Hand	Seat Depth	Weight
KM3520F16	16"	26"				20"	16"	37.5"	16"	29 lbs.
KM3520F18	18"	28"				20"	16"	37.5"	16"	29 lbs.
802SOMA-16	16"	25 1/4"	12 1/2"	34"	43 1/2"	19 1/2"	17"	37 1/4"	16 1/2"	29 lbs.
802SOMA-18	18"	26 1/2"	12 1/2"	34"	43 1/2"	19 1/2"	17"	37 1/4"	16 1/2"	29 lbs.
802F16-DY	16"	25"	12"	33"	42 1/2"	20"	18"	37"	16"	30 lbs.
802F18-DY	18"	26"	12"	33"	41 1/2"	20"	18"	37"	16"	31 lbs.
KM-8020F16	16"	24"	14"	30 1/2"	41"	20"	16"	37 1/2"	16"	30 lbs.
KM-8020F18	18"	26"	14"	30 1/2"	41"	20"	16"	37 1/2"	16"	30 lbs.
KM-8520-16	16"	24"	14"	30 1/2"	41"	20"	16"	37 1/2"	16"	30 lbs.
KM-8520-18	18"	26"	14"	30 1/2"	41"	20"	16"	37 1/2"	16"	30 lbs.

Q Quick Release Wheels F Flip-Back Armrests

Pediatric

MODEL No.	Seatwidth	Overallwidth	Foldedwidth	Length	Overall len.	Seat-Floor	Backrest	Floor-Hand	Seat Depth	Weight
BU1	13 1/2"	20 1/2"	16"	26 1/2"	34"	18"	22"	43 1/2"	15"	20 lbs.
KM-7520Q	13 1/2"	22 1/2"	12 1/2"	29"	33"	18"	15 1/2"	38"	13 1/2"	28 lbs.
KM-7501-TP	13 1/2"	20"	8"	22"	27 1/2"	15 1/2"	14"	40"	12"	20 1/2 lbs

