



Ultra Lightweight Adjustable Wheelchair

- Medicare code K0005 pending
- Padded seat and back sling
- Ergonomic hand rims
- Frame tubing 1 inch provides strength and rigidity
- Swing In & Away footrest
- "Tube-in Center" foot-plate, assures better side leg support
- Pouch for carrying small items attached to upholstery located behind backrest
- Backrest angle and back height adjustable
- Front seat and back height adjustable
- Multiple seat width and depth available
- 24"x1" Quick Release Rear spoke wheels w/polyurethane, high tread, flat free wheels std
- 8"x1" flat free front casters std
- Frame Color: Diamond Black
- Weighs 28.5 lbs (w/out footrest)
- Weight Capacity: 286 lbs
- Warranty: Frame Limited Lifetime





	seat	seat	armrest	seat	back	overall	overall	folded	overall	weight w/o	weight
	width	depth	height	height	height	height	open width	width	length	riggings	capacity
Model No.	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(lb.)	(lb.)
KM-8522Q16	16	16, 18	8-11	15-21	17/20/15-21	34-38	24	15	43	28.5	286
KM-8522Q18	18	16, 18	8-11	15-21	17/20/15-21	34-38	26	15	43	28.5	286
KM-8522Q20	20	16, 18	8-11	15-21	17/20/15-21	34-38	28	15	43	28.5	286