TRANZIT GO

FOLDABLE POWER WHEELCHAIR



PRODUCT FEATURES

- Weighs only 42 pounds (without battery)
- · Folding and lightweight frame
- · Padded seat and padded backrest
- · Flip up footplate
- Storage pouch under seat
- · Padded flip back armrest
- 7" x 1.5" solid polyurethane front caster
- 8" x 2" solid polyurethane rear caster
- Folded dimensions 24 x 11 x 31 inches
- Lithium ion battery (optional second battery)
- Removable joystick (mounted left or right side)
- Weight with one battery is 46 pounds (51 pounds w/ optional second battery)
- Turning radius 31.5"
- · Maximum obstruction height 1.20 inch
- · Maximum gradient slope 12 degrees
- · Weight capacity: 264 pounds

SPECIFICATIONS

Weight Capacity: 264 lbs

Front Casters: 7" x 1.5"

Rear Wheel: 8" x 2"

Battery: 24V/10AH Li-Ion

(2nd battery optional)

Max Speed: 3.75 mph (6 km/h)

Charging Period: 6-8 hours

Driving Range: 12 miles

Turning Radius: 31.5"

Overall Length: 39 inches

Overall Width: 24 inches

Overall Height: 36 inches

Folded Overall Width: 11 inches

Seat Width: 16 inches

Seat Depth: 16 inches

Back Height: 14 inches

Floor to Seat Height: 18.5 inches



