

XO-55

Horizon Standing Wheelchair Manual Propel

PRODUCT FEATURES



- 51 lbs total weight with wheels
- 6061 T-6 Aircraft Grade Aluminium Frame
- Pure white frame color
- Transparent lightweight composite side panel
- High strength struts for manual standing
- Rear strut for stability and increased nimbleness
- Easy to use standing handle allows for manual stand in as little as 10 seconds
- Quick release detachable rear wheels and caster forks
- Breathable nylon sling/upholstery seat and back
- 24" x 1" high profile pneumatic tires
- 6" x 1 1/4" poly foam front flat free casters
- Padded flip-back armrest and flip up footplates
- Chest strap, pelvic seatbelt, and knee straps are standard
- Scissor style wheel
- Lifetime warranty on frames

SPECIFICATIONS

HCPCS: N/A Seat Width: 16"/18" Overall Height: 32"
Seat Depth: 18/19/20 inch
Overall Width: 29"/31" Seat Height: 21"
Overall Length: 40" Back Height: 11"
Armrest Height: 8" Weight Capacity: 250 lbs



19255 San Jose Ave • City of Industry, CA 91792 •
(626) 581-2235 • www.karmanhealthcare.com