

XO-505

POWER STANDING WHEELCHAIR



PRODUCT FEATURES

- 6061 T-6 Aircraft-grade Aluminum Frame
- Full control LCD display
- Full leg rest support cushion
- Extendable arm/hand support bars on both sides
- H-style chest harness
- Dual front LED headlights (luxury model)
- Rear dual headlights w/ flashing feature (luxury model)
- New powered hydraulic system
- Full size frontal food tray (optional)
- Range per Charge: ~25 Miles
- Easy to use joystick for drive & stand
- USB support for mobile devices (optional)
- Power reclining feature
- Power tilting feature
- Mid wheel drive
- Full sit to stand power function

SPECIFICATIONS

HCPCS: N/A

Overall Height: 56.5 inch w/ headrest (46.6 inch without headrest)

Overall Width: 28 inch

Overall Length: 45 inch.

Armrest Height: 8.5 inch

Seat Width: 18 inch

Seat Depth: 18/19/20 inch

Seat Height: 23 inch

Back Height: 25 inch

Weight Capacity: 250 lbs

