



KARMAN HEALTHCARE
 19255 San Jose Ave.
 City Of Industry, CA 91748
 Tel: 626-581-2235 Fax: 626-581-2335



QR Reader Smartphone Info



XO-55 SERIES ORDER FORM

Submitting for QUOTE ORDER

Mark For: _____
 Date: _____
 Dealer Acct#: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: (____) _____ FAX: (____) _____
 Confirmation Email: _____

Purchase Order #: _____
 Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: _____

The HCPCS CODES are not intended to be, nor should be considered billing or legal advice. They are provided based on code verification through PDAC or our interpretations of billing. The use of HCPCS CODES does not ensure coverage or payment. For coverage information, please verify the appropriate payer's coverage policy

XO-55 Series:

XO-55 Lightweight Manual Propel Power Standing Wheelchair std features: flip back padded armrest, chest strap, pelvic seat belt, knee wedge support, 24"x1" flat free rear wheels, 6"x1 1/4" flat free front casters, two piece flip up footplate, calf support, anti tippers, 55 lbs w/battery, weight capacity 250 lbs.

BASE FRAME				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	XO-55N		16" Manual Propel Manual Propel and Standing Wheelchair, White	\$9,450.00
<input type="checkbox"/>	XO-55N-L		16" w/ extended seat depth Manual Propel and Standing Wheelchair, White	\$10,140.00
<input type="checkbox"/>	XO-55		18" Manual Propel and Standing Wheelchair	\$9,450.00
<input type="checkbox"/>	XO-55-L		18" w/ extended seat depth Manual Propel and Standing Wheelchair, White	\$10,140.00
SEAT DEPTH				
<input type="checkbox"/>	SDXO55-18		18" front to rear of seat sling	Std
<input type="checkbox"/>	SDXO55-19		19" front to rear of seat sling	Std
<input type="checkbox"/>	SDXO55-20		20" front to rear of seat sling	Std
<input type="checkbox"/>	SD-XO55-21		21" front to rear of seat sling, for "L" only	Std
<input type="checkbox"/>	SD-XO55-22		22" front to rear of seat sling, for "L" only	Std
SEAT CUSHIONS				
<input type="checkbox"/>	CU-ERGO-18	E2601	18x17 seat memory foam 2" cushion	\$176.00
<input type="checkbox"/>	CU-ERGO-16	E2601	16x17 seat memory foam 2" cushion	\$176.00
<input type="checkbox"/>	CU-FT-18	E2601	18x16 seat foam cushion	\$72.00
<input type="checkbox"/>	CU-FT-16	E2601	16x16 seat foam cushion	\$72.00
<input type="checkbox"/>	CU-GFR-18	E2601	18x16 seat gel cushion	\$98.00

<input type="checkbox"/>	CU-GFR-16	E2601	18x16 seat gel cushion	\$98.00
--------------------------	-----------	-------	------------------------	---------

BACK CUSHIONS				
<input checked="" type="checkbox"/>	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	CU-ERGO-18-back	E2611	18" wide memory foam 2" back cushion	\$116.00
<input type="checkbox"/>	CU-ERGO-16-back	E2611	16" wide memory foam 2" back cushion	\$116.00
<input type="checkbox"/>	BKF-1817	E2611	18" wide back foam cushion	\$98.00
<input type="checkbox"/>	BKF-1617	E2611	16" wide back foam cushion	\$98.00

REAR WHEELS				
<input checked="" type="checkbox"/>	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	RW-XO55		24x1 flat free polyurethane quick release wheel assembly	Std

KENDA TIRE COLOR/SPINERGY SPOKE COLOR				
<input type="checkbox"/>	KENDARED		Red Kenda tire	Std

OPTIONS & ACCESSORIES				
<input type="checkbox"/>	AT-XO55		Anti Tippers-pair	Std
<input type="checkbox"/>	LS16	K0038	16" width leg strap-calf support	\$19.00
<input type="checkbox"/>	LS18	K0038	18" width leg strap-calf support	\$19.00
<input type="checkbox"/>	CP2	E1399	Universal carry pouch-small	\$18.00
<input type="checkbox"/>	TD-15		Tie Down Kit Foldable 15ft/4.5m	\$39.99